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« 3D HIP PLANNER SOFTWARE »

www.oneortho-medical.com

ACCESS TO :

FRENCH..... 3

ENGLISH..... 69

Fabriqué par

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www.oneortho-medical.com

Table des matières

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www.oneortho-medical.com

Table of content :

ENGLISH.....	69
1. GOOD TO KNOW	71
2. GENERAL.....	75
3. TERMS OF USE.....	79
4. LIFETIME OF THE DEVICE	80
5. PLANNING ACCESS.....	81
6. GENERAL DISPLAY.....	82
7. SURGICAL PREFERENCES.....	85
8. VALIDATING A PLANIFICATION	87
8.1. BACKGROUND ON THE COURSE OF A CASE	87
8.2. DISPLAY/AUTOMATIC CALCULATIONS	93
8.2.1. Main window of the « Case” step.....	93
8.2.2. Main window of the “Cup” stage.....	98
8.2.3. Main window of the "Stem" step.....	100
8.2.4. Main window of the "Summary" step.....	105
8.3. PLANNING IMPLANTS AND ASSOCIATED SECTIONS	109
8.3.1. Planification settings: General.....	109
8.3.2. Dimension et positionnement of the cup.....	113
8.3.3. Stem sizing and positioning	115
8.3.4. Synthesis and adjustment of the femoral neck.....	118
8.4. DISPLAY TOOLS	119
8.5. PLANNING VALIDATION.....	123
9. APPENDICES	126
APPENDIX 1 : ANATOMICAL ENTITIES	126
10. MAINTENANCE.....	132
11. CYBERSECURITY.....	132
12. INCIDENTS / BUGS.....	132
13. CONTACT	132
14. FABRICANT	133
15. DEVICE INFORMATION AND ASSOCIATED SYMBOLS	133


GOOD TO KNOW

1. GOOD TO KNOW

USE OF THE SOFTWARE :

Purpose of the software
Logout


The purpose of the software is to schedule hip surgeries, especially total hip replacements.
It is important to log out after scheduling a case(s).
To do this, click on your name at the top right, then click "Logout".


GOOD TO KNOW

WHAT YOU NEED TO KNOW ABOUT THE "CUP" TAB

I want to change the tilt value


I want to change the anteversion value


GOOD TO KNOW

I want to change the size of the cup


I can't view value information

GOOD TO KNOW

WHAT YOU NEED TO KNOW ABOUT THE "STEM" TAB

The template displayed is not a representation of the patient's case being planned. It is an explanatory schematic model that allows you to understand what the values correspond to.

2. GENERAL

NAME / TRADE NAME OF THE DEVICE:

The name of the device is « 3D HIP PLANNER ».

INTENDED USE:

The 3D HIP PLANNER is intended for hip surgery planning in a pre-operative context.

GENERAL DESCRIPTION:

The 3D HIP PLANNER enables 3D visualization of hip and proximal femur anatomy; it allows the planning of the positioning, size, and model of a total hip prosthesis, in preparation for hip surgery. It includes a 3D library with several implant models, in their different sizes, of total hip prostheses. Once the planning has been performed and validated by the surgeon user, the software provides a planning report.

Using 3D items reconstructed from the DICOM images provided by the surgeon, this software displays the patient's hip and proximal femur bone structure. It displays the hip prosthesis within the digitally reconstructed bone structure. The surgeon user can manipulate parameters such as the range and size of the prosthesis. The default position of the prosthesis is set by One Ortho in such a way as to reduce the case's processing time for the surgeon user and is not a planning proposal.

The software's visual interface is available in dark and light versions.

SOFTWARE VERSION:

The current version in this manual is version 2.4.

PRECAUTIONS / WARNINGS:

- Before using the device, go through the User Manual: all steps and tools are described.
- Always follow good practices related to the surgery, see section **“Risk control measures external to the software”** stated below.
- Never transfer access rights to a third party, login and password are strictly personal.
- Make sure to understand the label. In case of doubt, consult the detailed information in the User Manual.
- Ensure that the version of the User Manual used is identical to the version identified on the software label.
- Always follow good practices related to the surgery, see section **“Risk control measures external to the software”** stated below.
- Always use a browser compliant with the manufacturer’s requirements and ensure that it is not obsolete in order to avoid security failures.
- Always follow good practices related to the surgery, see section **“Risk control measures external to the software”** stated below.
- Follow minimum requirements concerning hardware (including screen definition) and IT networks characteristics as defined in the User Manual to ensure correct display of device’s interface as to avoid inconvenience and security issues.
- To avoid issues with the planning output data, make sure the patient’s data is correct.
- Always follow good practices related to the surgery, see section **“Risk control measures external to the software”** stated below.
- Before validating the planning, check that the 3D visual of the synthesis page complies with the expected result. An unsuitable planning may lead to longer operating time.
- No temporary saving is possible, progress will be lost if the software is closed without having fully completed and validated the planning. Make sure to validate the case to save the planning and keep an archived planning report. After validation, the planning report is available on the platform.
- Please let us know if you experiment any bug (see §contact)
- At the user's request, the User Manual can be supplied in paper form at no additional cost to the user within 7 days of the request.

The use of the software allows a 3D planning carried out pre-operatively by the surgeon user¹. The data presented by the software are not a planning suggestion and the final choice of planning is made by the surgeon user².

The surgeon is required to have the patient's consent to process their personal data on the planning software.

It is forbidden to use multiple web browsers or multiple tabs of a browser in parallel to access multiple pages of the application.

GENERAL

INSTALLATION SPECIFICATIONS:

Operating Systems:

- Windows
- Mac OS.

Browser :

- Google Chrome

The browser used must be compatible with WebGL technology:

<https://caniuse.com/?search=webgl>

for 3D rendering quality. the only browser to use, which is verified and guaranteed by the manufacturer, is "Chrome".

There is no minimum as to the quality of the internet connection used but exchanges and transfers can be long with a low-speed connection. We recommend a connection with a throughput greater than 100Mbps.

The site is totally responsive, that is to say that it automatically adapts to the screen size of your hardware support with a minimum resolution of 1200x675.

UNITS USED


All lengths are indicated in millimetres (e.g. offset, length, etc.).

All angles are in degree (e.g. anteversion, tilt, etc.).

ACCESSIBILITY TO THE DOCUMENTATION:


The documentation is accessible from the platform at any time in the menu at the top right by clicking on the symbol next to the user's profile.

Then a menu appears offering to download the documents available according to your needs and the title above.

GENERAL

This User Manual is also accessible at any time within the application by clicking on the icon  .

The label above represents the validated version for this entire document.

TERMS OF USE & LIFETIME OF THE DEVICE

3. TERMS OF USE

INTENDED USE :

The 3D HIP PLANNER is intended for hip surgery planning in a pre-operative context. It will allow the determination of the adequate prostheses ranges and sizes to be implanted for a total hip arthroplasty.

INDICATIONS

3D HIP PLANNER is indicated when a total hip prosthesis is indicated.

INTENDED USERS

The intended users of the medical device are orthopaedic surgeons, with knowledge of anatomy, biomechanics and reconstructive surgery of the musculoskeletal system, and associated surgical techniques. The surgeon is responsible for planning and approving the plan with the medical device.

TARGET POPULATION

The 3D HIP PLANNING SOFTWARE is intended for use on males or females, having reached skeletal bone maturation and requiring total hip surgery as stated in the indications. Pregnant and breast-feeding women are excluded from the target population.

PERFORMANCES CLAIMED

No clinical performance is claimed, the claimed performances are the following:

- Assessment of needed stem size: a difference in size between the implanted prosthesis and the size planned when using the device of 1 or less, in 90% of cases.
- Assessment of needed acetabular cup size: a difference in size between the implanted prosthesis and the size planned when using the device of 1 or less, in 90% of cases.

BENEFITS

No clinical benefits are directly claimed for this software.

Clinical benefits are indirect and result of the total hip prosthesis.

Intended benefits are:

- Patient functional recovery regarding hip mobility,
- Pain reduction

TERMS OF USE & LIFETIME OF THE DEVICE

CONTRAINDICATIONS

The planning software must not be used:

- For joints other than the hip joint.
- For planning other prostheses ranges than those referenced in the library.
- When hip prosthesis is contraindicated.

UNDESIRABLE SIDE EFFECTS

No adverse side effects have been recorded in association with the device.

4. LIFETIME OF THE DEVICE


The software's lifetime is two years.

During this period, One Ortho guarantees the software's performance and undertakes to carry out updates to enable its optimal use.


PLANNING ACCESS

5. PLANNING ACCESS


From the dashboard, the user selects the **tab referring to the 3D HIP**

A new page appears:

They then have access to all the information relating to ongoing interventions. This is also where the creation of a new tracker is done. To do this, the user


must click on

GENERAL DISPLAY

6. GENERAL DISPLAY

In this section, we will describe all the information indicated and the possible interactions in this insert common to all the pages of the application. Depending on the need, the user will be able to navigate from the "CASE" page to the "SUMMARY" page through the "CUP" and "STEM" pages.


Regardless of the page displayed, the top banner of the application looks like this:

This sequence of pages to be gone through via the arrow  on the top ribbon allows each step to perform either a control operation, a planning operation, or both.

At the top left of the application, the following surgery information is shown:

Reminder of the
lateralization case
L for Ledt knee
R for right knee


Surgery information : Patient's and
date of birth, date of surgery,
Lateralization reminder and
patient's unique identification
number

Tools specific to the realization of specific actions are also available (measurement only on SUMMARY page):

GENERAL DISPLAY

Enabling and disabling the point to point leasurement tool

Modification of the colorimetry

Enabling and disabling comment

Pressing the Colorimetry button

will display two sliders that allow you to adjust the contrast and brightness:


Slider to control the contrast

Slider to control the brightness

Validation tools are available (only active on the SUMMARY page):

Final validation of the planning


Access to the planning report

Access to the planning report is active only after the planification is validated

Finally, tools for display management, access to explanatory documentation and disconnection from the platform are available at any time during the planning:

GENERAL DISPLAY

It is possible to zoom in and out with the mouse wheel on the imagery displayed throughout the application.

Depending on your equipment mouse, trackpad, touchpad the action may be different to achieve the goal.

For example, on PC with touchpad it is necessary to make a "pinch open",

For the same action, on Mac with touchpad, you have to perform a "swipe with two fingers".


SURGICAL PREFERENCES

7. SURGICAL PREFERENCES

The user has the possibility to specify his preferred parameters for the realization of his planning.

To do this, he must access the preferences screen setting via the button .


The interface of this page is as follows:

The user will then be able to configure the following elements:

The appearance of the planning theme:

Light theme

Dark theme

SURGICAL PREFERENCES

- The appearance of the offset and length information insert:

Standard

i OFFSET : 51 / -6.7 / 53
 LENGTH : 801 / +3.9 / 801

Simplified

i OFFSET : -6.7 +
 LENGTH : +3.9

- The display of the axial view on the default stem page:

Default View

Default masking

- The default first display page:

Default "CAS" page

ONEORThO CASE CUP STEM SUMMARY
 NAME: (blank) RIGHT DATE OF SURGERY: 2018-10-20 CASE NUMBER: 10100-RHT DATE OF BIRTH: 10-01-1975 OPERATING ROOM: 0001
 TOTAL LENGTH: 801 / 801
 FEMUR ACETABULUM
 DIAMETER: 46mm / 50mm
 ANTEVERSION: 13° / 16°
 INCLINATION: 53°
 BONY OFFSET: 51mm
 BONY CDD: 118°
 441 / 442
 377 / 375
 1-ANTERIOR SURFACE
 CE

Default "CUP" page

ONEORThO CASE CUP STEM SUMMARY
 NAME: (blank) RIGHT DATE OF SURGERY: 2018-10-20 CASE NUMBER: 10100-RHT DATE OF BIRTH: 10-01-1975 OPERATING ROOM: 0001
 SIZE: 50 / 50
 INCLINATION: 53° / 53°
 ANTEVERSION: 16° / 16°
 ONEORThO Ceramic cementless
 LENGTH: 801 / +3.9 / 801
 OFFSET: 51 / -6.7 / 53

- Automatic or custom configuration of anteversion and inclination values
- The default information for Cup and Stem implants:
 - The default range
 - The default type

To validate any changes made and return to the planning, the user must click on the button



VALIDATING A PLANIFICATION

8. VALIDATING A PLANIFICATION

8.1. BACKGROUND ON THE COURSE OF A CASE

As mentioned above, the tool was developed with the succession of 4 pages (screens) that must be gone through sequentially from the initial screen of the "Case" to the final page of "Synthesis" passing through the planning of the "Cup" and the planning of the "Stem".

Écran ① : 

On a typical computer screen, the initial window named « ① - Cases » is as follows for a surgery planned on the **right side** indicated by the pictogram  at the top left.

This first screen is a step of possible access to the user manual of the device and a step of control by the surgeon of the administrative and anatomical data of the case.

VALIDATING A PLANIFICATION

The information on this first page will remain accessible on the following screens by clicking on the double arrow under the indication "OPERATED SIDE".

The screenshot displays the OneOrtho medical software interface. At the top, there is a navigation bar with the OneOrtho logo, a 'CASE' tab, and a 'CUP' tab. Below this, there are two identical case detail cards. Each card features a large 'R' in a dark square, indicating the right side of the patient. The case details include: NAME : DEMO-Right, DATE OF SURGERY : 20-08-2025, CASE NUMBER : DEMO-Right, DATE OF BIRTH : 01/01/1975, and OPERATED SIDE : Right. A green double arrow icon is highlighted under the 'OPERATED SIDE' field in both cards. Below the case details, there is a section for 'FEMUR ACETABULUM' measurements. On the left, there is a 'Frontal' X-ray image of a hip joint with a green circle and crosshairs indicating the measurement area. Below the X-ray, there is an 'Axial' view. The measurements are listed in a table:

FEMUR ACETABULUM	
DIAMETER	46 mm 50 mm
ANTEVERSION	13° 16°
INCLINATION	53°
BONY OFFSET	51 mm
BONY CCD	118°


Écran ② :  CASE CUP STEM SUMMARY >

The second window named "② - Cup" is as follows for a surgery planned on the **right side** indicated by the pictogram  at the top left.

VALIDATING A PLANIFICATION

This second screen is a planning step of the prosthetic cup in its bone environment. This step makes it possible to fix the final prosthetic joint center at the center of the planned cup.

Planning is carried out in spatial positioning, dimensional and implant typology. All changes made to the prosthetic cup in the left interactive insert are updated in real time in the 3D dynamic insert on the right.


The offset (negative value) corresponds to the medialization of the cup.

VALIDATING A PLANIFICATION


Écran 3 :  CASE CUP **STEM** SUMMARY >

The third window named "3 - Stem" is as follows for a surgery planned on the **right side** indicated by the pictogram **R** at the top left.


By clicking on the button  an optional axial section scanner view can be displayed or hidden, the rest of the page is unchanged. This display is also configurable by default in surgical preferences.

VALIDATING A PLANIFICATION


This third screen corresponds to the planning of the prosthetic femoral stem in its bone environment.

Planning is carried out in spatial positioning, dimensional and implant typology. All changes to the prosthetic stem in the left interactive insert are updated in real time in the 3D dynamic insert on the right.

On a typical computer screen, the last window named "4 - Synthesis" is as follows for a surgery planned on the **right side** indicated by the pictogram **R** at the top left.

VALIDATING A PLANIFICATION

his last screen is a final control step by the surgeon of the administrative data of the case, the planned prosthetic data, the digital data of the osteotomy and the functional data (length and offset) which can still be modified by the choice of the length of the neck.


Other modifications are always possible by returning to the "② - Cup" step or to the "③ - Stem" step. In the case of a return to the "② - Cup" step, the sequential passage by the "③ - Stem" step will be mandatory before returning to the "④ - Synthesis" screen.

DISPLAY / AUTOMATIC CALCULATIONS

8.2. DISPLAY/AUTOMATIC CALCULATIONS

8.2.1. Main window of the « Case” step

This window is composed of two large parts.

he **Static insert for displaying the case's anatomical bone data** that displays the case's anatomical bone data:


- The side, by the turning the operated side **R** on the left window in bright - see the case above. For a left side, the display will be as follows:

DISPLAY / AUTOMATIC CALCULATIONS

The calculation of the lengths (in mm) of the bone segments of the two lower limbs to assess the inequality of length between the two limbs and the origin of this inequality by segment with the following:

The femoral length is calculated between the center of the femoral head and the top of the tibial spines. The grayed value opposite indicates the same measurement on the opposite, non-operated side.


The length of the lower limb is calculated between the radiological U point of the hip and the center of the malleoli. The grayed value opposite indicates the same measurement on the opposite, non-operated side.

The tibial length is calculated between the top of the tibial spines and the center of the malleoli. The grayed value opposite indicates the same measurement on the opposite, non-operated side.


DISPLAY / AUTOMATIC CALCULATIONS

- The diameter (in mm) of the femoral head, as that of the acetabulum, is calculated automatically by superimposing a sphere as congruent as possible to the bone element concerned. The diameter of the sphere thus positioned gives the diameter of the bone element.


- The anteversion (in degree) of the femur is calculated between the axis of the femoral neck and the posterior bicondylar plane. For the acetabular cup, it is the inclination of the cotyloid plane measured in a horizontal plane with respect to the sagittal plane.

DISPLAY / AUTOMATIC CALCULATIONS


- The inclination (in degree) of the cotyloid plane is measured using the frontal plane relative to the horizontal plane.


- The femoral bone offset (in mm) is measured between the center point of the femoral head and the top point of the great trochanter and projected into the front plane.

DISPLAY / AUTOMATIC CALCULATIONS

- The CCD angle (in degree) is the angle between the diaphyseal axis and the neck axis measured in a front plane

The second Dynamic **3D visualization insert area** displays, on the side to be operated, the reconstruction of the bones in 3D with implants placed. Before any user intervention, the position of the sections and the size of the implants are those defined by the ONE ORTHO team at the "**Pre-planning**" stage (see § Accessibility to Documentation to consult the "Workflow User Manual")). The icons for the visualization tools

DISPLAY / AUTOMATIC CALCULATIONS


are described in paragraph "8.4 Display tools" and each icon has a tooltip on mouse-over.

In addition to these specific tools, the usual functions that can be used in 3D scenes with your pointer are active:

- With a mouse, with the left click to rotate the objects in the 3D scene, the right click held to translate the object into the 3D scene and finally the wheel to zoom +/-.
- With a trackpad, with the click held with one finger to rotate the objects in the 3D scene, the click held with two fingers to translate the object into the 3D scene and finally drag it with two fingers to zoom +/-.


8.2.2. Main window of the "Cup" stage

This window is composed of three main parts.


The **Static Case Anatomical Data Display Box** displays the anatomical case bone data in white and the prosthetic data in dynamics of changes in color:

DISPLAY / AUTOMATIC CALCULATIONS


- The inclination (in degree) of the cotyloid plane in white measured between this plane and the anatomical horizontal plane.
- The anteversion (in degree) of the cotyloid cup in color measured between the equatorial plane of the prosthetic cup and the anatomical sagittal plane.
- The inclination (in degree) of the cotyloid cup in color measured between the equatorial plane of the prosthetic cup and the anatomical horizontal plane.

The second **Dynamic 3D visualization insert area** displays the hip in 3D, with the cup being placed on the side to be operated. The 3D display is dynamic of the changes made by the user.

The icons for the visualization  tools are described in paragraph "8.4 Display tools" and each icon has a tooltip on mouse-over.

This button allows the red display of the milled surfaces between the prosthetic cup and the unmilled bone acetabulum.


In addition to these specific tools, the usual functions that can be used in 3D scenes with your pointer are active:

- With a mouse, with the left click to rotate the objects in the 3D scene, the right click held to translate the object into the 3D scene and finally the wheel to zoom +/-.

DISPLAY / AUTOMATIC CALCULATIONS

- With a trackpad, with the click held with one finger to rotate the objects in the 3D scene, the click held with two fingers to translate the object into the 3D scene and finally drag it with two fingers to zoom +/-.

The last **Interactive inserting area for planning, positioning and dimensioning of the acetabular cup implant** is described in the paragraph "8.3.2 Sizing and positioning of the cup".

In this area and in all the cup planning scanner views two points are shown:


- A fuschia dot, the same color as the cup, that represents the prosthetic joint center.
- A yellow dot that represents the center of the bony femoral head.
- A white dot that represents the center of the acetabulum.
- A cyan dot that represent the articular center of the stem.

In the initial position, the fuschia and white points merge.

8.2.3. Main window of the "Stem" step


This window is composed of three main parts.

DISPLAY / AUTOMATIC CALCULATIONS


Static display insert for planned bone and prosthetic anatomical data

The **Static Insert area for displaying bone anatomical data** of the case displays the location data of the prosthetic stem in its bone environment and the extraarticular values of length and offset.

DISPLAY / AUTOMATIC CALCULATIONS

- In blue, the distance measured along the diaphyseal axis between the head center and the top of the great trochanter.
- In purple, the distance measured vertically between the shoulder of the prosthesis and the bottom of the saddle of the bone neck.
- In green, the distance measured in the axis of the neck between the cutting plane of the neck and the top of the bony head.
- In pink, the distance measured perpendicular to the cut of the neck between the cutting plane of the neck and the upper base of the small trochanter.

- In blue, the values of length and bone offset on the operated side.
- In green, the prosthetic correction values to the bone length and offset values on the initial operated side.
- In white, the length and offset values on the opposite side.


he second **3D Dynamic Display Insert area** displays the femur in 3D, with the stem placed on the side to be operated. The 3D display is dynamic of the changes made by the user. The icons

for the visualization  tools are described in paragraph "8.4 Display tools" and each icon has a tooltip on mouseover.

This button allows the red display of the grated surfaces between the prosthetic stem and the ungrated bone femoral barrel.

DISPLAY / AUTOMATIC CALCULATIONS


In addition to these specific tools, the usual functions that can be used in 3D scenes with your pointer are active:

- With a mouse, with the left click to rotate the objects in the 3D scene, the right click held to translate the object into the 3D scene and finally the wheel to zoom +/-.
- With a trackpad, with the click held with one finger to rotate the objects in the 3D scene, the click held with two fingers to translate the object into the 3D scene and finally drag it with two fingers to zoom +/-.

The last Interactive **Femoral Implant Planning, Moving and Sizing Insert** area is described in "8.3.3 Stem Sizing and Positioning".

DISPLAY / AUTOMATIC CALCULATIONS


In this area and in all stem planning scanner views are shown four points:

- A cyan blue dot of the same color as the stem that represents the articular center of the stem.
- A white dot that represents the center of the bony femoral head.
- A fuschia dot that represents the center of the prosthetic cup planned in the previous step.
- A green dot that materializes the upper base of the small trochanter

DISPLAY / AUTOMATIC CALCULATIONS

8.2.4. Main window of the "Summary" step

This window is composed of two large parts.

The **Static Case Anatomical Data Display Box** that displays the anatomical case data box:


- The side, by the turning the operated side on the left window in bright - see the case above. For a right side, the display will be as follows:

DISPLAY / AUTOMATIC CALCULATIONS

- The calculation of the lengths (in mm) of the bone segments of the two lower limbs to assess the inequality of length between the two limbs and the origin of this inequality by segment with:

The femoral length is calculated between the center of the femoral head and the top of the tibial spines. The grayed value opposite indicates the same measurement on the opposite, non-operated side.


The length of the lower limb is calculated between the radiological U point of the hip and the center of the malleoli. The grayed value opposite indicates the same measurement on the opposite, non-operated side.


top of the tibial spines and the center of the malleoli. The grayed value opposite indicates the same measurement on the opposite, non-operated side.

DISPLAY / AUTOMATIC CALCULATIONS


- The diameter (in mm), the name of the range and the version of the planned cup:

- The size, range name and planned stem version:

- Bone positioning measurements of the stem:

- In blue, the distance measured along the diaphyseal axis between the head center and the top of the great trochanter.
- In purple, the distance measured vertically between the shoulder of the prosthesis and the bottom of the saddle of the bone neck.
- In green, the distance measured in the axis of the neck between the cutting plane of the neck and the top of the bony head.
- In pink, the distance measured perpendicular to the cut of the neck between the cutting plane of the neck and the upper base of the small trochanter.

DISPLAY / AUTOMATIC CALCULATIONS

The offset and length measurements resulting from the planning:

- In blue, the values of length and bone offset on the operated side.
- In green, the prosthetic correction values to the bone length and offset values on the initial operated side.
- In white, the length and offset values on the opposite side.

The second area, **Dynamic 3D visualization insert**, displays the 3D reconstruction of the bones, with the implants positioned on the side to be operated on according to the surgeon's plan.

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he icons for the visualization  tools are described in paragraph "8.4 Display tools" and each icon has a tooltip on mouseover. In addition to these specific tools, the usual functions that can be used in 3D scenes with your pointer are active:

- With a mouse, with the left click to rotate the objects in the 3D scene, the right click held to translate the object into the 3D scene and finally the wheel to zoom +/-.
- With a trackpad, with the click held with one finger to rotate the objects in the 3D scene, the click held with two fingers to translate the object into the 3D scene and finally drag it with two fingers to zoom +/-.


The last Interactive **Insert area for choosing the length of the final neck** is described in paragraph "8.3.4 Synthesis and adjustment of the femoral neck".

PLANNING IMPLANTS

8.3. PLANNING IMPLANTS AND ASSOCIATED SECTIONS

8.3.1. Planification settings: General

After loading the case, the first screen is displayed:

This display is non-modifiable and presents the pre-planning carried out by the ONE ORTHO teams.

On the left, the limb lengths of the two lower limbs to assess the inequality of initial length overall or by segment.

In the center, the bone anatomical values of the patient.


On the right, a 3D view centered on the pathological bone situation.

By clicking , the default planification screen of the acetabular cup is displayed and the "Cup" button changes to color to indicate that this step is in progress.

PLANNING IMPLANTS


This display is interactive and presents the pre-planning carried out by the ONE ORTHO teams. On the left the 2D displays on the scanner images (at the top of the front, at the bottom in axial section in the center in sagittal view) with dedicated pads per view for the movements of the cup, in the center and at the top the pad for choosing versions and sizes and on the right the real-time 3D display of the modifications (detail in 8.3.2 below).

By clicking , the default stem planification screen appears and the "Stem" button changes to color to indicate that this step is in progress.

PLANNING IMPLANTS


This display is interactive and presents the pre-planning carried out by the ONE ORTHO teams at the opening. On the left, the 2D displays on the scanner images (at the top of the profile, at the bottom in axial section in the center in frontal view) with dedicated pads for stem movements, in the center and top the choice of versions and sizes, in the center and bottom the calculated values of stem tracking and on the right the real-time 3D display of the modifications made (detail in 8.3.3 below).

By clicking on , the summary screen is displayed and the "Summary" button changes to color to indicate that this step is in progress.

PLANNING IMPLANTS


This display is interactive and presents the coincidentally assembly of the femoral stem attached to the femur and the center of the cup itself attached to the pelvis, this corresponds to the intraoperative reduction of the stem in the cup.

In the right and left views, the values compared between the pre-operation member, the post-operation member, and the opposite member are displayed. The only variable that can be modified here is the head diameter and neck length of the head used (detail in 8.3.4 below).

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
The user can :

- Change the size of the cup that has been pre-planned. To do this, in the central inset one must click on "+" or "-" and the modified value is displayed in the center as well as the 3D cup of the chosen size on all the views.

- Choose the cup range to use by clicking on it.


- Adjust the position of the cup in antero-posterior and mid-lateral with the two pads present in the frontal and axial views. Outer arrows of the pads while watching the movement of the cup on the scanner views to assess the desired position.

Frontal


et Axial pad

- Adjust the position of the cup in tilt and ante/retro version with the 4 central circular arrows of the pads while watching the movement of the cup on the scanner views to evaluate the desired position.

In the 2D views of the cup superimposed on the scanner cuts, it is possible to scroll through the scanner sections with the green slicer on the right to refine the desired position. The cross-section position displayed in one view is marked in other views by a white line. This marking makes it possible to control the position of the scanner cut considered. At the opening all the cuts pass through the articular center of the cup. The center of the cup is represented by the fuschia dot, the yellow dot materializes the center of the bony head

PLANNING IMPLANTS



The three logos at the bottom allow:

Return to the initial display.



To display the cup in 2D section (below left) or in 3D (below right).

PLANNING IMPLANTS

Display the bone in 3D or in CT slices.

8.3.3. Stem sizing and positioning

The user can :

- Change the size of the stem that has been pre-planned.
- To do this, in the central inset you must click on "+" or "-" and the modified value is displayed in the center as well as the 3D stem of the chosen size on all the views.
- For the type click on "STD" for the standard version, and on "VARA" for the varized version.
- For bony head select the diameter and neck length according to the expectations on length and offset.

PLANNING IMPLANTS

the movement of the stem on the scanner views to evaluate the desired position.

Adjust the position of the stem in antero-posterior (In the sagittal view) and mid-lateral (In the frontal view) with the 4 arrows located in the pads at the bottom of the concerned views (white boxes) while watching the movement of the stem on the scanner views to evaluate the desired position.

- Adjust the position of the stem in varus / valgus (In the frontal view), in Ante-Retro version (In the axial view) and in inclination (In the sagittal view) with the 6 circular arrows located in the pads at the bottom of the relevant views (white boxes) while watching

In 2D views of the stem with scanner sections superimposed, it's possible to scroll through the scanner sections with the green slicer to fine-tune the desired position. The scanner cut selected in one view is indicated in the other views by a white line. This allows you to check the position of the selected scan section.

The three logos at the bottom allow, as for the cup:


Return to the initial display.


Display the stem in 2D section (below left) or 3D section (below right).

PLANNING IMPLANTS

Display the bone in 3D or in CT slices.



PLANNING IMPLANTS

8.3.4. Synthesis and adjustment of the femoral neck

The user can :

- Change the neck length to adjust the limb length and/or final offset. To do this, in the insert of the implants, you must select the desired neck length (Short, Medium or Long) by clicking

on the buttons







- Take a measurement (in mm) between two points in the 3D scene by double-clicking to insert a first point and double-clicking a second time to insert a second point, a measurement in mm taken between the two points will be displayed. There are several measurements that can be taken.

To clear the display of measurements, click on the button this action will delete all the displayed measurements.

DISPLAY TOOLS

8.4. DISPLAY TOOLS

	 Reset view Back to the initial display, tooltip
	Show, hide or transparently display the right femoral bone by  Show / Hide / Transparent successive clicks, tooltip: In transparent mode the stem must be in Display mode.
	Show, hide or show transparency the pelvis by successive clicks,  Show / Hide / Transparent tooltip: In transparent mode the cup must be in Display mode.
	Show, hide or show in transparency the left femoral bone by  Show / Hide / Transparent successive clicks, tooltip: In transparent mode the stem must be in Display mode.
	Show, hide or show in transparency the bone cotyloid plane by  Acetabular plane successive clicks, tooltip:
	Show, hide or show in transparency the anatomical landmarks by  Anatomical landmarks successive clicks, tooltip:

DISPLAY TOOLS

	Display the bone, by successive clicks, in turn in CT-scan or 3D 3D View section, tooltip: _____
	Display the cup, by successive clicks, alternately in 2D or 3D section, Implant View tooltip:
	Display the stem, by successive clicks, alternately in 2D or 3D section, Implant View tooltip::
	Display the head, by successive clicks, alternately in 2D or 3D section, Implant View tooltip::
	Display the areas of intersection between the implant and the bone, tooltip: « Afficher/Masquer surface fraisée (en rouge) » on the Cup page et « Afficher/Masquer surface râpée (en rouge) » on the Stem page. In this mode the pelvis or femur is in transparent mode and the implants must be in display mode to evaluate the milled or grated area
	Add an additional view
	Hides a view
	Allows you to scroll through the scanner images in the view where the slicer is present.

DISPLAY TOOLS

Frontal Axial Sagittal	Indicates the nature of the anatomical view
M	Indicates anatomical positioning M for Medial
L	Indicates anatomical positioning L for lateral
S	Indicates anatomical positioning S for Superior
I	Indicates anatomical positioning I for Lower
A	Indicates anatomical positioning A for Anterior
P	Indicates anatomical positioning P for Posterior
441 / 442 	Provides contextual information about the associated value
	By clicking gives information about the associated label.




DISPLAY TOOLS


	Clicking minimizes the display
	Clicking expands the display

PLANNING VALIDATION


8.5. PLANNING VALIDATION

When all the planning elements meet the user's expectations, he must position the 3D view on the right-hand side of the overview page according to the view he wishes to find on the second page of the planning report. This is indicated by a tooltip on the 3D view, which disappears as soon as the 3D view is clicked.

Select , an information insert appears indicating that the planification validation is final and that a planification report will be generated with the final view defined above.

PLANNING VALIDATION

To validate the planification, select **Confirm**, returns to the summary screen and always allows you to adjust the parameters.

Cancel

returns to the summary

After validation, the custom planning report appears and can be downloaded

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

PLANNING VALIDATION

This two-page document contains all the information of the case with all the anatomical and prosthetic information, on the second page is inserted the view defined by the user.

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

9. APPENDICES



APPENDIX 1 : ANATOMICAL ENTITIES



THE PUBIC SYMPHYSIS	
Définition	The area between the coxal bones. It is a fibrous cartilage tissue.
Taking the point	We choose the lowest area at the level of the pubic arch.
Good to know	This point will be defined as the origin of the 3D scene
	

ISCHIAL TUBEROSITIS	
Definition	These are the 2 lowest points in the frontal plane
Taking points	Select the extreme part of the ischium (right then left);
What you need to know	These are the points of contact with the horizontal plane if the pelvis were to be placed on this plane.
	

THE CENTER OF THE ACETABULA	
Definition	The acetabulum is a joint cavity of the iliac bone, located on each side of the pelvis, in which the femoral head (hemispherical end of the femur) is articulated to form the hip joint.
Taking the point	Select the inner section of the cotyloid horn

PLANNING VALIDATION


	 Then automatic creation of a sphere congruent to this area, determination of the center of this sphere as the center of the acetabulum and measurement of the diameter as the diameter of the bone acetabulum.
	

THE BOTTOM OF THE PSOAS NOTCH	
Definition	Lowest point of the curvature when looking at the acetabulum in the sagittal plane
	

PLANNING VALIDATION


THE ACETALUTE PLAN

Definition	We use three points to determine the plane of the acetabulum that will be used for tilt and anteversion angle measurements.
Taking points	Select the edge of the aceta to create an average plane of the selected surface.

THE SACROILIAC LOWER POINTS


Definition	The sacroiliac joints define the articular areas between the sacrum and the two iliac wings.
Taking the point	Choose the lowest area at the level of the pubic arc.

PLANNING VALIDATION


THE ANTERIOR POINTS OF THE RADIOLOGICAL "U"

Definition	<i>They are located in the inferomedial direction in the acetabulum, just above the obturator hole</i>
Taking the point	<i>The pelvis is cut through a frontal plane that passes through the centers of two acetabulums. Keep the posterior part and place the points of the radiological U.</i>

THE PROXIMAL END OF THE GREAT TROCHANTER

Definition	<i>The proximal end of the greater trochanter is the highest point of the largest protuberance on the upper part of the femur.</i>
------------	--


THE UPPER BASE OF THE PETIT TROCHANTER

PLANNING VALIDATION

Definition	<i>The femur has a second, smaller protuberance called the small trochanter, the upper part of the small trochanter is in front view of where a change in curvature is observed.</i>


THE POINT OF CURVATURE OF THE NECK	
Definition	<i>It is the point of curvature, in the plane of the face of the ridge that joins the great trochanter at the femoral head.</i>

THE DIGITAL DIMPLE	
Definition	<i>The digital dimple is located in the hollow between the inner surface of the greater trochanter and the femoral neck.</i>

POSTERIOR CONDYLE POINTS	
Definition	<i>These are the most posterior points in the sagittal plane</i>
 Points Postérieurs condyle	

THE CENTER OF THE FEMORAL HEAD	
Definition	<i>The femoral head can be modeled by a sphere. It is the center of this sphere</i>
Point taking	<i>Select the outer surface of the femoral head. Then automatic creation of a sphere congruent to this surface, determination of the center of this sphere as the center of the head and measurement of the diameter as the diameter of the bony femoral head.</i>

PLANNING VALIDATION

TOP OF THE FEMORAL HEAD	
Definition	The furthest point in the direction of the pass
	

MAINTENANCE, CYBERSECURITY & CONTACT

10. MAINTENANCE

An annual update is carried out for safety maintenance, automatically and completely transparently to the user. This maintenance includes resolving bugs and updating dependencies such as the framework, function libraries and SOUPs. It aims to prevent the emergence of new bugs and security vulnerabilities.

Depending on their criticality, some bugs can be resolved without delay.

11. CYBERSECURITY

Our approach to managing cybersecurity is based on risk management at each stage of the software lifecycle, integrating cybersecurity concepts right the design stage. Security measures are applied to user authentication to the web application, which is managed by a login and password that gives only one user access to his or her own cases, thus restricting access to patient data to authorized persons only. We therefore ask you to be careful not to divulge your password, and to log out when you stop using the application.

All cybersecurity risks are integrated and controlled in the medical device's risk management file.

12. INCIDENTS / BUGS

Despite the rigorous development of this software, any serious incident in connection with the device must be notified to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

A list of known bugs is available at the following link: [Link residuals bugs](#)

If you detect a bug that is not on this list, please send the information to the ONE ORTHO service.

13. CONTACT

Contact :

For any questions concerning the use of the platform, please contact ONE ORTHO at the following email address: contact@oneortho-medical.com or by phone at 04 26 78 76 74.

The ONE ORTHO team is available to answer you from Monday to Friday from 9am to 6pm.

MAINTENANCE, CYBERSECURITY & CONTACT

14. FABRICANT

ONE ORTHO
Parc INOPOLIS,
206 Route de Vourles
69230 Saint Genis Laval
FRANCE
contact@oneortho-medical.com

Class Im medical device, manufactured by the company ONE ORTHO,
complies with the requirements of EU Regulation 93/42/CEE.

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15. DEVICE INFORMATION AND ASSOCIATED SYMBOLS

Light theme:

 Manufacturer	 Date of manufacture	 Unique device identifier	 Medical device
 Consult the User Manual			

Dark theme:

 Manufacturer	 Date of manufacture	 Unique device identifier	 Medical device
 Consult the User Manual			

General Procedure case

Copyright Information

This document and what it describes are protected by the law of March 11, 1957 on literary and artistic property, modified by the law of July 3, 1985, as well as by copyright laws and international conventions.

Without the written permission of ONE ORTHO no part of this manual, including the products and software described herein, may be reproduced, transcribed, stored in a database system, or translated into any language, in any form or by any means, except for documentation kept by the purchaser for backup purposes. Any person failing to comply with these provisions shall be guilty of counterfeiting and liable to the penalties provided by law.

ONE ORTHO provides this manual "as is" without warranty of any kind, express or implied, or conditions of merchantability or fitness for a particular purpose. In no event shall ONE ORTHO, its directors, officers, employees, or agents be liable for any indirect, special, incidental or consequential damages even if ONE ORTHO has been advised of the possibility of such damages resulting from any defect or error in this manual or product.

Product and company names appearing in this manual may or may not be registered trademarks or subject to copyrights of their respective companies and are used only for purposes of identification or explanation, and for the sole benefit of the owners, without intent of infringement.