MANUEL UTILISATEUR / USER  
MANUAL /  
MANUALE UTENTE / MANUAL  
DE USUARIO

« 2D HIP PLANNER SOFTWARE »

[www.oneortho-medical.com](http://www.oneortho-medical.com)

# ACCESS TO :

|                      |           |
|----------------------|-----------|
| <b>FRENCH.....</b>   | <b>3</b>  |
| <b>ENGLISH.....</b>  | <b>30</b> |
| <b>ITALIAN .....</b> | <b>57</b> |
| <b>SPANISH.....</b>  | <b>84</b> |

**FRENCH**

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**« 2D HIP PLANNER SOFTWARE »**

[www.oneortho-medical.com](http://www.oneortho-medical.com)

DOC 0035 H  
17 octobre 2025

Table des matières

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## COPYRIGHT INFORMATION

This document and the information contained herein are protected by the French law of 11 March 1957 on literary and artistic property, as amended by the law of 3 July 1985, as well as by copyright laws and international agreements.

No part of this manual, including the products and software described in it, may be reproduced, transcribed, stored in a database system or translated into any other language, in any form or by any means, apart from documentation retained by the buyer as a backup, without the written permission of OneOrtho Medical. Anyone breaching these rules will be guilty of copyright infringement and will face criminal penalties.

OneOrtho Medical provides this manual "as is" without any express or implied warranty of any kind including merchantability and fitness for a particular purpose. Under no circumstances shall OneOrtho Medical, its directors, managers, employees or agents be liable for indirect, special, accidental or consequential damage even if OneOrtho Medical was informed of the possibility of such damage resulting from faults or errors in the manual or product.


The names of the products and companies that appear in this manual may or may not be trademarks or subject to copyright for their respective companies, and are used only for identification or explanatory purposes and for the sole benefit of the owners, with no intended breach.

The specifications and information contained in this manual are provided for information only and are subject to change without prior notice and shall not be interpreted as a commitment from OneOrtho Medical. OneOrtho Medical is in no way liable for any errors or inaccuracies in this manual, including the products and software described herein.

The software is a decision-aid system designed for professionals having completed the appropriate medical training. It must in no way be used as the only basis for making clinical decisions during the diagnosis, care or treatment of a patient. The plausibility of the information obtained via this software must be systematically clinically checked before being used in the treatment of patients. Any use of the programme's medical information which is not included in the original concept or intended usage is not recommended and will be considered as improper use of the software.

A list of known bugs is available via the following link: [Residuals bugs](#)

CE  
0459

This device complies with the requirements of directive 93/42/EEC relating to medical devices.  
CE mark acquisition year: 2017

© 2017 OneOrtho Medical. All rights reserved.

# GENERAL

## 1. GOOD TO KNOW

For an optimal experience, ensure that your system is configured to display all platform information correctly.

### MAC USERS:

On macOS, you must enable the permanent display of scroll bars in order to view all implant proposals correctly.

To do this:

1. Open the **Apple menu** > **System Preferences**.
2. Select **Appearance**.
3. In the **Show scroll bars** section, choose **Always**.

## 2. GENERAL

### DEVICE NAME / TRADE NAME:

The name of the device range is "2D HIP PLANNER SOFTWARE".

### DESTINATION:

This product is designed for preoperative planning of hip surgery.

### VERSION OF THE SOFTWARE:

The version presented in this manual is 1.5.

### INSTALLATION SPECIFICATIONS :

#### Operating systems :

- Windows,
- Mac OS.

#### Browsers :

- Google Chrome
- Safari

The zoom level must be 100% for optimal use of the platform.

#### Screen definition:

The platform is fully responsive, meaning that it adapts to the size of the screen. On smaller formats such as smartphones, this adaptation can modify the layout of the user screen. It is recommended to use a minimum of 1200px in 16:9 to display all the information on the application.

#### Internet connection:

An internet connection is required to use the application. A minimum connection speed of 10 Mbps is recommended.

#### Hardware:

# GENERAL

The hardware required for the client workstation for a comfortable user experience is:

- Processor: The processor must be equivalent to an Intel Core i5 (sixth generation or later)
- Graphics card: The graphics card must have a minimum video memory of 2048 MB.

# GENERAL

## ACCESSIBILITY OF DOCUMENTATION:

This User Manual can be accessed at any time on the application by clicking on the icon

OneOrtho

medicall

R

LAST NAME :  
FIRST NAME :  
INTERVENTION DATE :

New

INFORMATIONS

LAST NAME \* :

FIRST NAME \* :

INTERVENTION DATE :

OPERATED SIDE \* :

PLANNING \* :

Choose a file

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Advanced planning:

Free planning:

OUTIL DESSIN

IMPLANTS

Drag and drop an image

DOC 0035 H  
October 17<sup>th</sup>, 2025

p 35/109

# INFORMATION

## 3. INFORMATION

To access 2D planning from the platform, click on the "2D planner" module in the dashboard.

The interface is arranged as follows:

The first step is to enter the patient information. Enter the last name and first name by clicking on

The top right banner consists of 4 buttons :

- The  button switches the planner to full screen mode.
- The  button allows the user to access the preferences.
- The  button allows you to save the user manual.
- The  button allows you to quit planning and return to the platform home page.

You can then choose the patient's side by clicking on Left.

# INFORMATION

PLANNING \* :

ADVANCED

FREE

Then select **"advanced"** or **"free"** planning by clicking on (further details in section §4. Advanced Planning or §5. Free Planning).

Choose a file

Then click on . A browser window opens in which you can choose X-ray images in jpg format and click on OK. The image is displayed again.

The information tab is displayed:

OneOrtho

R

LAST NAME : Test  
FIRST NAME : Demo  
INTERVENTION DATE : 2025-04-25

New

INFORMATIONS

LAST NAME \* :

FIRST NAME \* :

25/04/2025

INTERVENTION DATE :

OPERATED SIDE \* :

ADVANCED

FREE

PLANNING \* :

Choose a file

radio001.jpg

DOES THIS IMAGE HAVE A SINGLE HIP ? \*

YES

NO

CALIBRATE IMAGE

Calibration Explanation:

Click on the center of the ball and adjust the size of the ball circle.


OUTIL DESSIN

IMPLANTS

20 cm

# INFORMATION

You can then use the mouse wheel to zoom in on the image, particularly the ball which will then be used to scale the image.

CALIBRATE IMAGE


Then scale the image by clicking on .  
A window opens. Enter the size of the ball in mm and click on OK.

Enter the size of the ball in mm and draw a circle around it.

25 mm

Then, using the mouse, draw a circle around the ball to scale the image. Start in the middle of the ball and hold the mouse button down to enlarge the circle until it corresponds with the ball.

**Important:** The contours of the circle must be inside the contours of the ball or the known size marker.  
***Make sure you draw the circle by dragging the mouse from left to right to obtain positive values.***

Selecting the grey dot on the edge of the circle allows you to change its radius while the grey dot in the centre allows you to move the whole circle.

You then need to indicate whether the X-ray shows one or two hips by clicking on:

DOES THIS IMAGE HAVE A SINGLE HIP ? \*

YES

NO

OK

Click on "Confirm " to confirm and continue planning.


# "ADVANCED" PLANNING

## 4. "ADVANCED" PLANNING

This section describes how to use 2D planning after drawing notable geometric elements like the femoral shaft axis (intramedullary axis) and the acetabular cup centre. These elements will then be used for the automatic positioning of implants.

- [Drawing Tools](#)


The "**Drawing Tools**" tab opens automatically if you selected "**Advanced Planning**". You can choose different drawing tools in this tab:

- The **TRAPEZOID** tool is used to draw a trapezium in the femur using the mouse. The trapezium bases represent two levels of femoral shaft section.


The middles of these bases are used to automatically define the anatomical axis of the femur.

## "ADVANCED" PLANNING


To create the trapezium on the selected side in the "Information" section, stop the mouse at each trapezium corner (for around 1 second) and then move on to the next dot while holding the mouse button down.

You can click on the grey dots to adjust the edges of the trapezium or in the middle to move the whole drawing.


- The **LESSER TROCHANTER** tool allows you to measure the offset between a patient's two hips. For example, the difference in height in the lesser trochanter between the right and left sides. This tool will be used if there are two hips on the uploaded image.

## "ADVANCED" PLANNING

- The **CIRCLE** tool allows you to draw a circle from the centre of the bony acetabulum. Selecting the dot on the edge of the circle allows you to change its radius while the grey dot in the centre allows you to move the whole circle.


The  icon appears once these functions have been used. This enables you to delete the line drawn using the relevant function.


You can then click on  to continue .

- [Implants](#)


## "ADVANCED" PLANNING

OK

Clicking on  opens the **"Implant"** tab automatically. The femoral stem and acetabular cup implants and templates appear on the image.


You can choose which type of femoral stem you want to plan by selecting it from the drop-down menu in the **"Femoral stem"** section, then the type of acetabular cup in the same way in the **"Acetabular cup"** section.


Users can choose the type of pole they want to plan by clicking on the drop-down menu under the heading **"Stem"**.


# "ADVANCED" PLANNING


Then its type of cupule by clicking on the drop-down menu under the heading "Cup".


You can:

- Change the opacity of the drawing tools from 0% to 100% using this slider  for adjustment.
- Change the colour of the stem using the  buttons under the title "Stem", the first three buttons respectively serving to quickly colour the selected stem in predefined colors, and the 4th button used to open a colour picker to choose a precise colour of the user's choice.
- Move the femoral stem along the anatomical axis defined by the trapezium by clicking on the diagonal arrows.

# "ADVANCED" PLANNING

- Change the colour of the cup using the Color:    Other:  buttons under the title " Cup ", the first three buttons respectively serving to quickly colour the selected cup in predefined colors, and the 4<sup>th</sup> button used to open a colour picker to choose a precise colour of the user's choice.

Two additional items of information are given and displayed on the X-ray itself: offset and height.

The offset indicates the horizontal gap between the centre of the femoral head at 0 on the femoral stem and the centre of the acetabular cup.

The height indicates the vertical gap between the centre of the femoral head at 0 on the femoral stem and the centre of the acetabular cup.


It is also possible to plan implants on the opposite side using the buttons

The button on the selected side appears with a green background to indicate the selection

# "ADVANCED" PLANNING

When the side opposite the planning side is selected, an implant appears in the center of the image. The stem and the cup are then manipulated separately, as for those used for planning.


For these implants, the pad movement buttons are deactivated and the free movement button is available


It is then possible to drag and drop the implant layer onto the image to place it correctly. When the stem movement tool is activated, the cup movement tool is automatically deactivated and vice versa. To avoid a handling error by forgetting to deactivate the displacement tool, it deactivates automatically shortly after a displacement.

It is also possible to change the size of the implants on the directional pads of the stem and the cup, as well as the rotation of the implants without impacting the planning made on the side of the operation:

## "ADVANCED" PLANNING


Paddle "Stem" on the opposite side

"Cup" pad on the opposite side


# "ADVANCED" PLANNING

Click on  to continue

Clicking on  opens a window to save the file you have created in PDF format, with the information relating to the patient, the original image and the plan you have just created.

The following window is displayed:

By clicking on , the case will be saved in the case archive module. By clicking on , the case will not be saved.

Saving the case consists of saving the planning report in the form of a PDF file. This saving does not allow the elements of the planned case to be retrieved in order to modify it later. For any modifications, the case will have to be recreated and processed again according to the desired modifications.

# "FREE" PLANNING

## 5. "FREE" PLANNING

Choosing **"Free"** planning takes you directly to the **"Implant"** tab. You can choose which type of femoral stem you want to plan by selecting it from the drop-down menu in the **"Femoral stem"** section, then the type of acetabular cup in the same way in the **"Acetabular cup"** section.

Options are displayed as follows:

The screenshot shows the **IMPLANTS** planning interface. It is divided into two main sections for selecting components.

**Top Section (Acetabular Cup):**

- SYMBOL:** A drop-down menu showing "Cementless dual mobility" with a value of 48.
- Color:** Three colored circles (green, yellow, grey) and an "Other:" option with a green circle.
- Plan cup on opposite side:** A toggle switch.
- Display cup on opposite side:** A toggle switch.
- Navigation:** "LEFT" and "RIGHT" buttons, and a central circular diagram with arrows and the number 48.

**Bottom Section (Femoral Stem):**

- SYMBOL:** A drop-down menu showing "Cementless without collar" with a value of 1.
- STD / VARA:** Two buttons, with "STD" currently selected.
- Color:** Three colored circles (green, yellow, grey) and an "Other:" option with a green circle.
- Plan stem on opposite side:** A toggle switch.
- Display stem on opposite side:** A toggle switch.
- Navigation:** "LEFT" and "RIGHT" buttons, and a central circular diagram with arrows and the number 1.


An **OK** button is located at the bottom center of the interface.

Users can choose the type of pole they want to plan by clicking on the drop-down menu under the heading "Stem".

This close-up shows the femoral stem selection interface. It includes a drop-down menu with the text **STEMDEMO** and "Cementless without collar", a value of 7, and two buttons labeled **STD** and **VARA**.


## "FREE" PLANNING

Then select your type of cup by clicking on the drop-down menu under the heading "Cup".

You can:


- Change the colour of the stem using the Color: buttons under the title "Stem", the first three buttons respectively serving to quickly colour the selected stem in predefined colors, and the 4<sup>th</sup> button used to open a colour picker to choose a precise colour of the user's choice.

Pad « Stem » in Free Planning

- Change the colour of the cup using the Color: buttons under the title "Cup", the first three buttons respectively serving to quickly colour the selected cup in predefined colors, and the 4<sup>th</sup> button used to open a colour picker to choose a precise colour of the user's choice.


## "FREE" PLANNING

Pad « Cup » in Free planning

It is then possible to drag and drop the implant template onto the image to position it correctly. When the stem displacement tool is activated, the cup displacement tool is automatically deactivated and vice versa. To avoid a handling error by forgetting to deactivate the displacement tool, it deactivates automatically shortly after a displacement

Before being moved, the implants appear in the center of the image


It is also possible to place implants on the opposite side using the buttons

## "FREE" PLANNING


When the side opposite the planning is selected, the selected implant appears in the center of the image. The stem and the cup are then manipulated separately as for those used for planning.

For free planning of implants on the opposite side, the pads remain identical and it is possible to freely move

the layers using the button, to change the size of the implants of the stem and the cup, as well as the rotation of the implants without impacting the planning done on the side of the operation.


Click on  Validate case


## "FREE" PLANNING


Clicking on  opens two windows. One to save in PDF format the file that the user has just created, with the patient information, the original image as well as the schedule they just made:

And another window is opening

By clicking on , the case will be saved in the case archive module. By clicking on , the case will not be saved.


The case can be saved by saving the planning report in the form of a PDF file. This backup does not allow the elements of the planned case to be retrieved for subsequent modification. For any modifications, the case will have to be recreated and processed again according to the desired modifications.

# VARIOUS OPTIONS


## 6. VARIOUS OPTIONS

Using the contrast scroll bar  , the user can change the contrast of the image. To do this, he must click on the position indicator on the scroll bar and move his mouse from right to left while continuing to press his mouse.

Using the brightness scroll bar  , the user can change the brightness of the image. To do this, he must click on the position indicator on the scroll bar and move his mouse from right to left while continuing to press his mouse.


## VARIOUS OPTIONS


## VARIOUS OPTIONS


The user can also use the tool that is used to display the value of an angle and measure to display the value of a length.

It is possible to have 3 angle values displayed at the same time. For each corner, a button appears to be able to delete it.

The user can also use the tool that serves to measure a length

It is possible to have 3 lengths displayed at the same time. For each length, a button appears to be able to delete it.

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The screenshot displays the "IMPLANTS" planning interface. It is divided into two main sections for selecting components. The top section is for the "Cotile" (cup) and the bottom section is for the "Stelo" (stem).

**Top Section (Cotile):**

- SYMBOL:** A dropdown menu showing "Cementless dual mobility" with a value of 48.
- Color:** Three circular color selection buttons (dark green, light green, yellow) and an "Other:" button with a dark green circle.
- Plan cup on opposite side:** A toggle switch.
- Display cup on opposite side:** A toggle switch.
- Navigation:** "LEFT" and "RIGHT" buttons flanking a central circular navigation pad with arrows and a central value of 48.

**Bottom Section (Stelo):**

- SYMBOL:** A dropdown menu showing "Cementless without collar" with a value of 1.
- STD / VARA:** Two buttons for selecting the stem type.
- Color:** Three circular color selection buttons (dark green, light green, yellow) and an "Other:" button with a light green circle.
- Plan stem on opposite side:** A toggle switch.
- Display stem on opposite side:** A toggle switch.
- Navigation:** "LEFT" and "RIGHT" buttons flanking a central circular navigation pad with arrows and a central value of 1.

An "OK" button is located at the bottom center of the interface.

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This is a close-up of the stem selection menu. It features a dropdown menu with the text "STEMDEMO" and "Cementless without collar", showing a value of 7. Below the dropdown are two buttons: "STD" (highlighted in green) and "VARA".

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
The screenshot displays the 'OneOrtho' application interface. At the top, there is a header bar with the 'OneOrtho' logo, a green button with a white 'R', and fields for 'LAST NAME:', 'FIRST NAME:', and 'INTERVENTION DATE:'. To the right of these fields are icons for 'New', a clock, and a power button. Below the header, the main interface is divided into two sections. The left section, titled 'INFORMATIONS', contains a form with the following fields: 'LAST NAME \*:', 'FIRST NAME \*:', 'INTERVENTION DATE:', and 'OPERATED SIDE \*:'. The 'OPERATED SIDE \*:' field has two buttons, 'R' and 'L', with 'R' selected. Below these fields is a section for 'PLANNING \*:' with two buttons, 'ADVANCED' and 'FREE', with 'ADVANCED' selected. There is also a 'Choose a file' button and a text label 'Aucune image sélectionnée'. The right section is a large white area with a dashed border and a central icon of a square with an upward arrow, labeled 'Drag and drop an image'. At the bottom of the 'INFORMATIONS' section, there are two buttons: 'OUTIL DESSIN' and 'IMPLANTS'.

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The screenshot shows the OneOrtho software interface. At the top, there's a header with the OneOrtho logo, a patient ID 'R', and fields for LAST NAME, FIRST NAME, and INTERVENTION DATE. Below the header, the interface is split into two main sections. On the left is a sidebar with a dark green header 'INFORMATIONS' and a light gray body containing various input fields and buttons. On the right is a large central workspace with a light gray background, a white rectangular area, and a 'Drag and drop an image' prompt with an upload icon.

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A dark green rectangular button with the word 'INFORMATIONS' in white capital letters.

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A light gray form with two input fields for 'LAST NAME \*' and 'FIRST NAME \*'. Below these is a section for 'OPERATED SIDE \*' with two buttons, 'R' (Right) and 'L' (Left), where 'R' is highlighted in green.

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A light gray form with a section for 'PLANNING \*' containing two buttons: 'ADVANCED' (highlighted in green) and 'FREE' (gray).

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Choose a file

The screenshot shows the OneOrtho medical software interface. The top header includes the OneOrtho logo, a patient ID 'R', and fields for 'LAST NAME: Test', 'FIRST NAME: Demo', and 'INTERVENTION DATE: 2025-04-25'. On the right, there are icons for 'New', a double arrow, and system settings. The main area is divided into two panels. The left panel, titled 'INFORMATIONS', contains several input fields and buttons: 'LAST NAME \*' (Test), 'FIRST NAME \*' (Demo), 'INTERVENTION DATE' (25/04/2025), 'OPERATED SIDE \*' (R/L), a 'PLANNING \*' section with 'ADVANCED' (selected) and 'FREE' buttons, a 'Choose a file' button next to 'radio001.jpg' (highlighted with a red box), a 'DOES THIS IMAGE HAVE A SINGLE HIP ? \*' section with 'YES' and 'NO' buttons, a 'CALIBRATE IMAGE' button, and a 'Calibration Explanation' section. The right panel displays a large X-ray image of a pelvis with a 20 cm scale bar on the right side.

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This is a close-up screenshot of the 'INFORMATIONS' tab from the software interface. It shows the 'Choose a file' button next to the filename 'radio001.jpg'. Below this is the question 'DOES THIS IMAGE HAVE A SINGLE HIP ? \*' with 'YES' and 'NO' buttons. A dashed line separates this from the 'CALIBRATE IMAGE' button. At the bottom, the 'Calibration Explanation' section is visible, with the text 'Click on the center of the ball and adjust the size of the ball circle.'

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The screenshot shows the 'IMPLANTS' configuration screen. It features two main sections for selecting implants. The top section is for 'Cementless dual mobility' (Symbol 48) and the bottom section is for 'Cementless without collar' (Symbol 1). Each section includes a 3D model of the implant, a 'SYMBOL' dropdown menu, color selection options (Color and Other), and toggle switches for 'Plan cup on opposite side' and 'Display cup on opposite side'. Below these are 'LEFT' and 'RIGHT' buttons, a central directional pad, and an 'OK' button at the bottom.

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This close-up shows the 'STEMDEMO' dropdown menu, which is currently set to '7'. Below the dropdown are two buttons: 'STD' (highlighted in green) and 'VARA'.

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