## NOTICE D'UTILISATION: KNEE CUSTOM-MADE GUIDES - (OneCustom Knee PSI)

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## INSTRUCTIONS FOR USE: KNEE CUSTOM-MADE GUIDES - (OneCustom Knee PSI)

Before using a medical device manufactured by One Ortho, the surgeon must take note of the following recommendations:

### 1. IMPORTANTE

Custom-made guides should only be used by the defined users, for the indication intended by the manufacturer and only for a specific patient in conjunction with our customers' implant range.

The following points are of decisive importance for successful use:

- Practice according to the current state of science and surgical art,
- Knowledge of the product's specific surgical technique.

The quality of the results obtained can only be guaranteed if all the instructions given in this manual are known and followed.

### 2. GENERAL INFORMATION

**Product trade name:** OneCustom Knee PSI

**Product description:** Single-use custom-made guides are associated with the PTG (Total Knee Prosthesis) of One Ortho customers. These guides are used to position the proximal tibial cut, the distal femoral cut and the patellar cut. These are three separate guides, one for the femur, one for the tibia and one for the patella, which can be made available together or independently of each other, depending on demand. These guides are said to be made-to-measure, as they are adapted to the shape of the patient's bone (femur, tibia or patella) and can only be used with a single patient.

**Medical use:** Custom-made guides are designed for a specific patient in orthopedic knee surgery, to prepare the bone for correct positioning of a total knee implant. This enables the implant to be inserted.

### 3. MANUFACTURING AND MATERIALS

ONE ORTHO's made-to-measure guides are manufactured by 3D laser sintering in PA12 polyamide. They are designed and manufactured in line with the latest state-of-the-art surgical knowledge and rigorous quality criteria, both in the choice of materials used and in the manufacture of the instruments.

Each instrument is identified by its marking and labeling, which includes the batch number (associated with the patient), the instrument identification and a size marking. Materials used in the manufacture of instruments comply with French, European or international standards setting out the characteristics of medical-grade materials.

### 4. PRODUCT INDICATIONS

Custom-made knee guides are only used in conjunction with an implant range offered in the Digital Planner Total software for knee arthroplasty, and are therefore indicated when the associated implant range is indicated, i.e. in the following cases:

- Degenerative osteoarthritis or post-traumatic osteoarthritis, in patients whose age, weight and activity level are compatible with a satisfactory long-term result.
- Rheumatoid arthritis.
- Revision of unicompartmental arthroplasty procedures, if bone stock and ligament stability permit.
- Failed osteotomy.

### 5. INTENDED USERS

This medical device is intended for orthopedic surgeons with knowledge of anatomy, biomechanics and reconstructive surgery of the musculoskeletal system, as well as the corresponding surgical techniques. No additional training is required for its use. It is sufficient to read the manual and the corresponding surgical technique carefully.

### 6. TARGET POPULATION

The target population, with the exception of pregnant or breast-feeding women, consists exclusively of bone-mature adult patients (men and women) with acute pain and/or joint disease resulting from a medical condition for knee arthroplasty.

### 7. CLAIMED PERFORMANCE

One Ortho does not claim any clinical performance for custom-made knee guides.

The technical performance claims are as follows:

- Custom-made knee guides maintain their properties after sterilization, during use and up to the expiration date.
- Custom-made knee guides are compatible with the corresponding total knee prosthesis.
- Custom-made knee guides are correctly positioned and sized for prosthesis implantation.

### 8. CLINICAL BENEFITS

One Ortho does not claim any clinical benefit for custom-made knee guides.

The clinical benefit claimed for the patient is indirect: it results from the placement of the implant, which in turn has the following benefits:

- Relieve patient pain,
- Restore mobility,
- Reduced operating time compared to conventional instrumentation.

### 9. CONTRAINDICATIONS

Contraindications to the use of custom knee guides are as follows:

- Use on joints other than the knee joint,
- Unsuitability for the manufacturer's intended use,
- Use with implants that are not part of the same One Ortho customer product range,
- Contraindications for One Ortho's range of total knee prostheses.

### 10. SIDE EFFECTS AND POSSIBLE COMPLICATIONS

Side effects related to the use of custom-made knee guides are listed below:

- Post-operative infection
- Inflammatory response or an allergic reaction
- Possible tissue injury.

### 11. LIFETIME

Custom-made guides have a lifetime of just one use, and are disposable.

### 12. Shelf life

Shelf life for custom-made knee guides is indicated on the label in YYYY-MM format (see "Shelf life" symbol): it corresponds to the last day of the month indicated.

### 13. WARNINGS / PRECAUTIONS FOR USE

### NEVER REUSE A CUSTOM-MADE GUIDE.

It is the surgeon's responsibility to familiarize himself/herself with the manufacturer's recommendations, including the accompanying documents provided for the information of healthcare professionals, in order to identify custom-made guides and use them in the conditions defined by One Ortho.

The surgeon must familiarize himself/herself with the instruction leaflet and the surgical technique before using the custom-made guides.

The surgeon and/or the patient must inform ONE ORTHO and the competent authority of the member state in which they are established of any serious incident involving custom-made knee guides.

Custom-made knee guides may only be used with the range of implants indicated on the planning report supplied with the guide.

Information for users / residual risks and associated recommendations are listed below:

Risk: "Patient allergy"	To avoid this residual risk, please follow this recommendation: "Make sure the patient is not allergic to polyamide".
Risk: "Instrument contamination" Risk: "Non-sterile product" Risk: "Poorly sterilized instrument"	To avoid this residual risk, please follow this recommendation: "Follow the information in the instructions concerning the method and conditions of sterilization".
Risk: "Labelling information is not available on the product" Risk: "The information on the label is not correct (UDI / Date of production/...) prior to production / use of the guide".	To avoid this residual risk, please follow this recommendation: "Check compliance of custom-made knee guides before use".
Risk: "Missing instrument"	To avoid this residual risk, please follow this recommendation: "Ensure availability of instruments before starting surgery".
Risk: "Loss of patient/instrument concordance" Risk: "Reversal between 2 patients".	To avoid this residual risk, please follow this recommendation: "Check the conformity of the instrument before using it".
Risk: "Conflict between the femoral cutting guide and the patella" Risk: "Impossibility of drilling the fixation holes through the guides" Risk: "Poor bone preparation" Risk: "Impossible recutting" Risk: "Impossibility of removing the guide" Risk: "Cutting thickness error" Risk: "Non-precise control of the bone cut".	To avoid this residual risk, please follow this recommendation: "Observe the instructions for use given on the label, in the package leaflet and in the operating instructions".
Risk: "Contamination through contact with a soiled instrument".	To avoid this residual risk, please follow this recommendation: "Follow the information given in the instructions for waste disposal".
Risk: "Re-sterilization of a guide that has fallen on the floor"	To avoid this residual risk, please follow this recommendation: "Discard any instrument that has fallen on the floor without attempting to clean or re-sterilize it".

#### 14. STORAGE AND HANDLING

- Instruments must be stored in a place protected from moisture and heat.
- The functional surfaces of the instruments must not be damaged, impacted, abraded, or otherwise deteriorated. Deformations of the instruments can lead to improper positioning of the prosthesis. It is therefore recommended that the instruments be handled with the utmost care.
- Custom-made guides must be stored in their original sealed packaging. Check the expiry date of the instruments on the packaging. Reject any instrument whose packaging is not intact.

#### 15. STERILIZATION

Custom-made knee guides are delivered clean, non-sterile and suitable for autoclave sterilization. All guide references are packaged in medical paper bags that are compatible with autoclave sterilization, which must be carried out by the user before use.

To achieve a Sterility Assurance Level (SAL) of 10<sup>-6</sup> One Ortho recommends :

Pre-vacuum steam autoclave sterilization.	18 minutes	134° C	Recommended for all instruments (France)
---	------------	--------	--

The re-sterilization of custom-made knee guides is only permitted in exceptional circumstances, and only before their first use, in the following cases:

- If the expiration date has passed,
- In case of failure during the first sterilization cycle.

This resterilization must be performed without a second cleaning step.

Any instrument that has fallen to the floor must be discarded. It must not be cleaned or resterilized under any circumstances.

**NOTE 1: concerning French instruction DGS/R13/2011/449 of December 1, 2011** on recommendations to be observed during invasive procedures in order to reduce the risks of transmission of non-conventional transmissible agents: the level of risk presented by the use of ancillary equipment from the One Ortho company justifies the choice of its treatment by the usual procedure of pre-disinfection, cleaning and sterilization in a pre-vacuum steam autoclave with a sterilization tray lasting 18 minutes at 134°C. The use of sterilization containers adapted to the type of sterilization chosen is the sole responsibility of the healthcare establishment, particularly for maintenance and control purposes.

**NOTE 2:** Any use of detergent and/or disinfectant for custom-made guides prior to use is the responsibility of the user and is not recommended by One Ortho.

#### 16. INSTRUCTIONS FOR USE

- Full metal instrumentation must be available during surgery, throughout the use of custom-made knee guides.
- Installation and positioning of instruments
  - Instruments are designed to adapt to the patient's anatomy. Surfaces in contact with the instruments must be kept clear to ensure correct positioning. Take the necessary time to position them on the patient. Ensure that no anatomical structures have been damaged during positioning.
- Drying
  - Ensure that the guide does not move when drilling the fixation holes
  - Use appropriate equipment to avoid excessive heat build-up. Do not use a drill bit larger than the diameter of the instruments.
- Cutting
  - Make sure the guide does not move during cutting
  - Use appropriate equipment to avoid excessive heat build-up. Do not use blades thicker than the slots in the instruments.

Refer to the Operating Technique supplied by the manufacturer for proper use of custom-made knee guides.

#### 17. Elimination

The professional waste producer is responsible for sorting infectious risk waste and waste that can be assimilated to household waste. Generally speaking, ordinary waste such as packaging can be treated as household waste and follows the household waste circuit. Objects and materials soiled by blood and/or body fluids are treated as Waste from Healthcare Activities with Infectious Risks.

#### 18. DEVICE INFORMATION AND SYMBOLS

	Manufacturer		Catalogue number		Batch code		Unique identifier device		Medical Device
	Expiration date		Keep out of direct sunlight		Moisture sensitive		Non-sterile		Patient initials
	Consult instructions for use		Do not use if packaging is damaged and consult instructions for use		Do not reuse		Prescribing physician		Surgery date
	Associated case number								

Any serious incident involving custom-made knee guides must be reported to the manufacturer and to the competent authority of the state in which the user and/or patient is established, stating the name, reference and batch of the product and the nature of the incident or malfunction.

**Manufacturer**  
ONE ORTHO  
Parc INOPOLIS – 206 Route de Vourles  
69230 SAINT GENIS LAVAL, FRANCE  
Phone: +33 (0)4 26 78 76 74  
[contact@oneortho-medical.com](mailto:contact@oneortho-medical.com)  
[www.oneortho-medical.com](http://www.oneortho-medical.com)

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