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Before using a medical device manufactured by One Ortho, the surgeon must take note of the following recommendations:

1. IMPORTANT INFORMATION

Custom-made orthognathic expansion devices must only be used by the defined users, for the indication determined by the manufacturer and only for one specific patient.

- The following points are key to guaranteeing that these medical devices are used successfully:
- Surgical practices aligned with current scientific knowledge and state-of-the-art surgical techniques,


Knowledge of the specific operating technique required for the product.

The quality of the results obtained can only be guaranteed if all of the instructions set out in this document are known, understood, and followed

2. GENERAL INFORMATION

Product trade name: Custom-made bone-supported surgical distractor

Product description: The orthognathic expansion device is a bone-supported (hard palate) custom-made medical device that allows the palate to be widened at the maxilla. It is composed of a positioning guide and a surgical distractor, the guide making it easier to position the surgical distractor.

Medical purpose: Orthognathic expansion devices are intended for a specific patient, to be used as part of orthognathic surgery on the palate of the maxilla (mid-palatal expansion) in order to correct transverse bone abnormalities.

3. MANUFACTURING AND MATERIALS

ONE ORTHO orthognathic expansion devices are manufactured by 3D printing in a titanium alloy (TA6V ELI). Their design and manufacture take the most recent knowledge of state-of-the-art surgical techniques and practices into account and meet exacting quality criteria, both in the selection of materials used and in manufacturing practices.

The materials used for the manufacturing of the orthognathic expansion devices comply with French, European and international standards stipulating the specifications for medical-grade materials.

4. PRODUCT INDICATIONS

Selective orthognathic surgery on the maxilla in cases of a transverse bone abnormality.

5. INTENDED USERS

Maxillofacial surgeon for the surgical aspect and patient for device activation.

6. TARGET POPULATION

Adults and pediatric population.

7. CONTRAINDICATIONS

The orthognathic expansion device is contraindicated in the following patients:

- Patients with a flat and/or scarred cleft palate.
- Patients with gum disease or periodontal disease.
- Patients with poor oral hygiene.
- Patients with a history of immunodeficiency, corticosteroid therapy, coagulation disorders, uncontrolled endocrine pathology, rheumatic disease, bone disease, diabetes or cirrhosis, or any other systemic or acute disease.
- Patients with osteomyelitis or active infection.
- Patients with a metal allergy or foreign body sensitivity
- Patients who have received radiotherapy to the head
- Patients with poor circulation and insufficient bone structure (insufficient bone quantity) or potential bone defects (insufficient bone quality) in the area where the transpalatal distractor will be fitted

8. WARNINGS/PRECAUTIONS FOR USE

- No adjustments are tolerated with orthognathic expansion devices (guide included).
- *The product is designed for single use only. The reconditioning of a used product (recycling) and its reuse in a patient are not authorized.*

9. STORAGE AND HANDLING

- Orthognathic expansion devices should be in a cool, dry place.
- Orthognathic expansion devices should not be subject to any damage, impact, abrasion or other deterioration. Do not use the medical device if it is damaged.

10. Cleaning and sterilization

Orthognathic expansion devices are supplied non-sterile and should be cleaned and sterilized prior to surgical use. Before cleaning, remove all of the original packaging. Prior to steam sterilization, package the product using a packaging system that has been approved for this type of sterilization. Follow the cleaning and sterilization instructions below

1. AUTOMATIC CLEANING:

Equipment: Washer/disinfector in accordance with the standard NF EN ISO 15883-1 and -2, detergent (neodisher® MediClean — concentration 0.5%), rinsing product (neodisher® MediKlar — concentration 0.05%), demineralized water. Please note: Preferably, a manual cleaning method with additional machine-washing should be used.

1. If the instruments are jointed, make sure that they are wide open or, if they're in several parts, take them apart. Arrange the openings and cavities so that water can drain out.

2. Start a cycle:

- Cold water prewash without any product,
- Wash for at least 15 minutes at a temperature below 50°C (45°C),
- Rinse with cold demineralized water.

3. When unloading, inspect instrument surfaces and cavities to check that all visible dirt has been completely removed.

Automatic disinfection::

Program a final rinse with thermal disinfection (27 mins. / demineralized water at 93°C), the time depends on the size and power of the device (NF EN ISO 15 883-1).

Drying::

This step avoids bacterial proliferation, removes the marks left by the washing and rinsing processes, and optimizes sterilization efficacy. - When using a washer/disinfector, do not exceed 120°C (dry using air at 110°C, then cool with air at 30°C).

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NOTE 1 : Any use of detergent and/or disinfectants other than the one recommended for orthognathic expansion devices is at the sole responsibility of the user

11. INSTRUCTIONS FOR USE

Before surgery

- The orthognathic expansion device must be cleaned and sterilized before use.
- The instruments and associated accessories must be prepared before any surgical procedure: the orthognathic expansion device is compatible with osteosynthesis screws with a diameter of 2 mm and a screwdriver from the ORTRAUTEK range (GLOBAL D), a practitioner activation key and a patient activation key.

During Surgery

Fitting the orthognathic expansion device (figure 1) :

The orthognathic expansion device is fitted during surgery, before or after the bone is cut for the distraction (which allows the surgical field and the approaches to be cleared)

1. The metal positioning guide is secured to the rod of the surgical distractor and is positioned on the palate
2. The positioning of the rod and the axis of each screw are precisely determined in advance

Figure 1

Figure 2

Fixing (figure 2)

The orthognathic expansion device is secured using trans mucosal fixation with osteosynthesis screws with a diameter of 2 mm, from the ORTRAUTEK range from GLOBAL D. During the planning phase, the spatial positioning of the screw axis allows for the length of each screw to be determined and anticipated.

The 3D design of the eyelets will guide the screwing. It is recommended that the screws are inserted and tightened progressively, alternating among the 4 screws as indicated.


Figure 3

Removal of the positioning guide (figure 3)

Before finalizing the screw fixings, the guide must be removed. The holding ligatures are cut so that the surgeon can remove the positioning guide .

Figure 4


End of screw fixing and activation (figure 4)

Finish gently screwing in the screws and gradually pull the rod to rest it on the mucosal membrane, without compressing it. An initial activation is carried out when the positioning guide is removed, at the end of the procedure, using the practitioner's activation key

Then, using the patient's activation key, the patient carries out the activation themselves, at a rate of one to two activations per day, according to the prescription issued by the surgeon and the orthodontist.

After surgery

- Activation by the patient is carried out using a patient activation key. Before use, this key must be disinfected using a biocide intended for human hygiene and which is compatible with titanium alloy. The activation is carried out with a quarter-turn (one quarter-turn corresponds to a spacing of 0.2 mm).
- To carry out an activation, follow the diagram below. Insert the key into the visible hole in the screw, gently push the key towards the bottom of the oral cavity, until you hear a "click" (follow the direction of the arrow marked on the device), then turn the screw. Gently release the key from the screw to prevent the screw from turning back
- Certain foods, such as sticky and chewy candy, chewing gum and crusty bread should be avoided for the safety of the patient and to ensure the proper functioning of the device. All other hard foods, nuts and stone fruit are also to be avoided.
- It is recommended to cut food into small pieces before putting them in the mouth and to avoid eating too much at once. You also must not touch the device or poke utensils, such as cutlery (fork, knife, etc.) into your mouth's structures or into the device.
- The removal of the device and any modification or corrective actions must be carried out by a healthcare professional.
- In the event of any abnormalities or if you have any doubts, consult your orthodontist or surgeon immediately.

12. DISPOSAL

The professional producing the waste is responsible for sorting potentially infectious waste and similar waste from domestic waste. In general, ordinary waste such as packaging is comparable to domestic waste and should be disposed of in the same way as domestic waste. Items and materials soiled by blood and/or bodily fluids should be treated as infectious clinical waste (classified as "DASRI" in France).

13. DEVICE INFORMATION AND ASSOCIATED SYMBOLS

	Manufacturer		Catalog reference		Batch number		Medical Device		Health care center or doctor
	Keep away from sunlight		Date of manufacture		Keep dry		Patient name		Consult instructions for use
	Do not use if the packaging is damaged and refer to the instructions for use		Caution		Do not reuse		Non-sterile		Patient Number
	Date of surgery								

Any serious incident occurring in connection with the expanders must be notified to the manufacturer and to the competent authorities in the country where the user and/or patient is registered, reporting the product name, reference and batch, as well as the nature of the incident or malfunction.

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