Fabriqué par / manufactured by / prodotto da / fabricado por

MANUEL UTILISATEUR
USER MANUAL
MANUALE UTENTE
MANUAL DE USUARIO

2D HIP PLANNING SOFTWARE

www.oneortho-medical.com

DOC 0103 H
Apr/20/2026

Page 1 / 169

Table of languages

FR **3**
EN **45**
IT **86**
ES **128**

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2D HIP PLANNING SOFTWARE

www.oneortho-medical.com

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OneOrtho
PATIENT INFO CALIBRATION DRAWING IMPLANTS CUTS

R ANONYMIZED ID
LAST NAME *
FIRST NAME *

PATIENT INFO

Fill fields.
* Required fields

Advanced planning:
Stem and cup layer planning attached to geometric construction tools for acetabular center and diaphyseal axis.

Free planning:
Stem and cup planning by free movement of layers overlaying the radiograph.

DEMO - 2D HIP PLANNING SOFTWARE
ONE ORTHO
Parc SOPHUS,
206 Route de Tourtes
83200 SAINT-ETIENNE-LE-MOULIN,
France
contact@oneortho-medical.com

UDI XXXXXXXXXXXXXXXXXXXXXXXX
AAAA-MM Read the PDF:
www.oneortho-medical.com/white

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Case number
Side

PATIENT INFORMATION

Anonymous code / LAST
NAME
First name
Case number

SURGERY INFORMATION

Surgery date
Surgeon
Planning date & time

2D PLANNING SUMMARY


	CUP	LINER	STEM	HEAD
DESCRIPTION / REFERENCE				
SIZE				

OFFSET LENGTH


2D planning image

IMPLANT RESIZING COEFFICIENT: BALL SIZE:

CORRECTION AND RESECTION DATA

ORIGINAL IMAGE

Manufacturer information

Manufactured by

USER MANUAL

2D HIP PLANNING SOFTWARE

www.oneortho-medical.com

Restriction of use:

This software is designed for professionals (Orthopaedic surgeons). It should in no way be used as the only basis for clinical decision-making for diagnosis, care or management of the patient. The plausibility of the information obtained through the software should be systematically clinically verified before it is used for the treatment of patients. Any application of the medical information of the program that does not fit into the original concept or intended use of the program is discouraged and will be considered an abuse of the software.

1.	GENERALITIES	48
1.1.	GOOD TO KNOW	48
1.2.	NAME / TRADE NAME OF THE DEVICE	48
1.3.	INTENDED USE	48
1.4.	GENERAL DESCRIPTION	48
1.5.	SOFTWARE VERSION.....	48
1.6.	PRECAUTIONS / WARNINGS	49
1.7.	INSTALLATION SPECIFICATIONS	52
1.8.	ENVIRONMENT OF USE	52
1.9.	ACCESSIBILITY OF DOCUMENTATION:.....	53
2.	TERMS OF USE	54
3.	LIFESPAN	56
4.	ACCESS TO PLANNING	57
5.	DISPLAY INFORMATION	59
6.	SURGICAL PREFERENCES	63
7.	PLANNING REALIZATION	65
7.1.	PATIENT INFORMATION	65
7.2.	CALIBRATION	68
7.3.	DRAWING	71
7.4.	IMPLANTS.....	73
7.5.	CUTS	78
7.6.	VALIDATION AND REPORTING.....	80
8.	MAINTENANCE	81
9.	CYBERSECURITY	81
10.	CONTACT	81
11.	MANUFACTURER	82
12.	INFORMATION CONCERNING THE DEVICE AND ASSOCIATED SYMBOLS	82
13.	ANNEX	84
13.1.	ANNEX 1 : PLANNING REPORT	84

Generalities

1. Generalities

1.1. GOOD TO KNOW

- PURPOSE OF THE SOFTWARE

The purpose of the software is to provide information to orthopaedic surgeons to help them in the planning of hip surgery, in particular total hip replacements.

- CONFIDENTIALITY

It is important to log out at the end of using the software. To do this, click on your name (top right) and click on "Log out".

1.2. NAME / TRADE NAME OF THE DEVICE

The device's name is "2D HIP PLANNING SOFTWARE".

1.3. INTENDED USE

2D HIP PLANNING SOFTWARE is intended to inform clinical management by providing 2D planning tools for preoperative hip surgery planning.

1.4. GENERAL DESCRIPTION

The "2D HIP PLANNING SOFTWARE" is a device that allows:

- Perform 2D planning (Superposition) of the implant model defined with the patient's bone structure in the form of X-rays in **.jpeg or .png format with minimal definition of 512px**
- Determine the implant range and define the most appropriate size for a given patient.

Access to the software is made through a secure web platform. For the realization of the 2D planning, the surgeon can load his X-ray on the associated module and once his planning is done, he has the possibility to download an automatically generated planning report in PDF format.

The application also allows you to save and archive planning data by a surgeon subject to patient authorization.

1.5. SOFTWARE VERSION

The active version of the software is version **2.3**.

***Note:** 2D Hip Planning software is designed to provide healthcare professionals with information in order to achieve planning hip surgeries. It is a standalone software, currently available on the market only in version 2.3. There is no other version of the software available to the users. The software can only assist with surgical planning if the intended prosthesis for implantation is available in the software's 2D scheme library.*

Generalities

1.6. PRECAUTIONS / WARNINGS

In spite of our great care in reducing all risks, some residual risks remain. The residual risks, their associated recommendations and patient consequences are listed below:

Residual risks	Recommendations	Patient consequences
Software understanding: The user does not understand how to position the prosthesis on the screen	Before using the device, go through the User Manual: all steps and tools are described. Always follow good practices related to the surgery, see section “Risk control measures external to the software” stated below.	Longer operating time.
Access to the software by someone who should not have access rights.	Never transfer access rights to a third party, login and password are strictly personal.	No patient consequence.
The design of the labeling format or content is not carried out according to regulatory requirements (symbols and regulatory information in the software)	Make sure to understand the label. In case of doubt, consult the detailed information in the User Manual.	No patient consequence
Wrong User Manual version available on the IFU website leading to software use error	Ensure that the version of the User Manual used is identical to the version identified on the software label. Always follow good practices related to the surgery, see section “Risk control measures external to the software” stated below.	Longer operating time.

Generalities

<p>The user uses an outdated or unverified browser.</p> <p>A security flaw due to using an outdated browser.</p>	<p>Always use a browser compliant with the manufacturer's requirements and ensure that it is not obsolete in order to avoid security failures.</p> <p>Always follow good practices related to the surgery, see section "Risk control measures external to the software" stated below.</p>	<p>Longer operating time.</p>
<p>The user does not understand what is the required recommended minimum configuration (including connection speed) to run the software.</p>	<p>Follow minimum requirements concerning hardware (including screen definition) and IT networks characteristics as defined in the User Manual to ensure correct display of device's interface as to avoid inconvenience and security issues.</p>	<p>No patient consequence. The medical staff must use other means to perform the plan of the surgery.</p>
<p>Incorrect patient data</p>	<p>To avoid issues with the planning output data, make sure the patient's data is correct.</p> <p>Always follow good practices related to the surgery, see section "Risk control measures external to the software" stated below.</p>	<p>Error in the data.</p> <p>Longer operating time.</p>
<p>The user does not know how to save his planning before to leave the session and loses their file.</p>	<p>After validation, the planning report is available on the platform.</p>	<p>No patient consequence. The medical staff must use other means to perform the plan of the surgery.</p>
<p>The user does not know how to contact the support in case of software bugs.</p>	<p>Please let us know if you experiment any bug (see §10 Contact).</p>	<p>No patient consequence.</p>

Generalities

The user needs to use the user manual in paper form.	At the user's request, the User Manual can be supplied in paper form at no additional cost to the user.	No patient consequence.
--	---	-------------------------

Risk control measures external to the software:

The use of 2D HIP PLANNING SOFTWARE **does not replace** the usual intraoperative control measures.

In order to avoid increasing risks of intraoperative and postoperative complications, the user is required to follow risk control measures external to the software: the user must follow the operative technique associated with the planned hip prosthesis provided by the prosthesis manufacturer. In particular the trial steps must be respected: the sizes are confirmed intraoperatively with anatomical landmarks, rasp testing, trial implants.

The use of the software allows 2D planning carried out preoperatively by the surgeon user². The data presented by the software is not a planning suggestion and the choice of planning is made by the surgeon user².

The surgeon is required to have the patient's consent to process their personal data on the planning software.

It is forbidden to use multiple web browsers or multiple tabs of a browser in parallel to access multiple pages of the application.

² Who holds a degree in orthopedic surgery, cf. intended user in section 2Terms of Use.

Generalities

1.7. INSTALLATION SPECIFICATIONS

Operating Systems:

- Windows,
- Mac OS.

Browser :

- Google Chrome

The browser used must be compatible with WebGL technology:
<https://caniuse.com/?search=webgl>

Screen definition:

It is recommended to use a minimum of 1200px in 16:9 to display all application information.

Internet connection:

An internet connection is required to use the application. A minimum connection speed of 10 Mbps is recommended.

Hardware:

The client workstation hardware required for a comfortable user experience is:

- **Processor:** The processor must be equivalent to an Intel Core i5 (Sixth generation or newer)
- **Graphics card:** The graphics card must have a minimum video memory of 2048 MB.

1.8. Environment of use

The intended environments of use are:

An at-home office


An office at the hospital or clinic

There are no particular characteristics required for the place of use. Only the installation specifications to be found in “ §1.7 INSTALLATION SPECIFICATIONS “ must be met in order to use the device safely.

Generalities

1.9. ACCESSIBILITY OF DOCUMENTATION:

This User Manual is accessible at any time on the application by clicking on the icon

The label above represents the validated version for this entire document.

Terms of Use

2. Terms of Use

Indications

2D HIP PLANNING SOFTWARE is indicated when a primary total hip prosthesis is indicated.

Total hip prosthesis may be indicated primarily in case of:

- primary or secondary osteoarthritis,
- femoral neck fracture,
- osteonecrosis or avascular necrosis,
- bone cancer.

Intended users

The intended users of the medical device are orthopaedic surgeons, with knowledge of anatomy, biomechanics and reconstructive surgery of the musculoskeletal system, and associated surgical techniques. The surgeon is responsible for planning and approving the plan with the medical device.

There is no specific software training required to use the device for the first time.

Target population

The 2D HIP PLANNING SOFTWARE is intended for use in patients requiring primary total hip arthroplasty. This population include both male and female adults of advanced age who have reached bone maturity (average age close to 60-70 years).

Pregnant and breast-feeding women are excluded from the target population.

Claimed performance

2D HIP PLANNING SOFTWARE gives information to clinical management by providing 2D planning tools for achieving accurate preoperative hip surgery planning:

- Planned and actually implanted prostheses sizes match to within one size in at least 75% of cases and within two sizes in at least 90% of cases.

Clinical Benefits

The device helps reduce both intraoperative and postoperative complications.

Contraindications

The planning software must not be used:

- For joints other than the hip joint and other models than those referenced in the library
- When hip prosthesis is contraindicated.
- On patients that have not reached skeletal bone maturation.

Terms of Use

Undesirable side effects

No adverse side effects related to the web application are listed.

Lifespan

3. Lifespan


This software has a lifespan of 2 years.

During this period, One Ortho guarantees the performance of the software and undertakes to make updates to allow its optimal use.


Access to planning

4. Access to planning

Click on the 2D HIP PLANNING module

A new page is displayed:


Then, he has access to all the information relating to the interventions in progress. This is also where the creation of a new planning is done. To do this, the user

must click on 

A new planning page is displayed :

Access to planning



Display Information

5. Display Information

In this section, we will describe all the information indicated and the possible interactions in this insert common to all the pages of the application. Depending on the need, the user can navigate from the "Patient Info" page to the "Cuts" page via the "Calibration", "Drawing" and "Implants" pages.

The top banner of the application looks like this, depending on the current page:


To move from one page to another, simply click on the navigation arrow  at the top right or click directly on the page navigation arrows:

-  ↔ The "PATIENT INFO" page is active.
-  ↔ The "CALIBRATION" page is active.

At the top left of the application, the following intervention information is shown:

Reminder of the lateralization of the case

L for a left hip
R for a right hip


Patient information: Anonymized code **OR** Surname/First Name of the patient must be entered

Display Information

Tools specific to the realization of specific actions are also available (only usable on the IMPLANTS and CUT pages):

Enabling and disabling the point-to-point measurement tool

Enabling and disabling the differential measurement tool

Activating and disabling the angle measurement tool


Enabling and disabling comment

Create the point-to-point measure: Click 2 points on the layer. The measurement will be displayed automatically

Delete the measure: Hover the mouse over the measurement drawing.  appears.


Click .




Create the angle measure: Click 3 points on the layer. The measurement will be displayed automatically

Delete the measure: Hover the mouse over the measurement drawing.  appears.


Click .

Display Information

	<p><u>Create the differential measure:</u> Click 2 points on the layer. Move the 3rd point created to display the differential measurement.</p> <p><u>Delete the measure:</u> Hover the mouse over the measurement drawing.  appears. Click .</p>	
	<p><u>Create an annotation:</u> Add your comment instead of "Enter your comment".</p> <p><u>Adding a new annotation:</u> Click  at the top right of the existing comment. Add your comment in the new window created.</p> <p><u>Move comment:</u> Click in the text box and, while pressing, move the comment.</p> <p><u>Delete an annotation:</u> Click  at the top right of the existing comment. Confirm deletion via the message that appears</p>	

Validation-specific tools are available:


Final validation is active only on the IMPLANTS and CUT pages.

If temporary backup is selected, the user will be redirected to intervention tracking. He can return to the current planning at any time by clicking on PLANNING in the "Planning" column of the desired case.

Finally, tools for display management, access to documentation and disconnection from the platform are available at any time during the planning process.:

Display Information

It is possible to zoom in and out with the wheel – mouse wheel on the imagery displayed throughout the application.

Depending on your equipment mouse, trackpad, touchpad the action may be different to achieve the goal.

For example, on PC with touchpad it is necessary to make a "pinch open",

For the same action, on Mac with touchpad, you have to perform a "swipe with two fingers".


Surgical Preferences

6. Surgical preferences

The user has the possibility to specify his preferred parameters for the realization of his planning.

To do this, he must access the screen setting his preferences via the button :

The interface of this page looks like this:


Surgical Preferences


The user will then be able to configure the following elements:

- The appearance of the planning theme:

Light Theme

Dark Theme

- The size of the calibration ball in mm (by default, it is 25.4 mm).
- The default planning mode: Advanced or free
- The default information for Cupule and Stem implants:
 - The default range
 - The default variant

To validate any changes made and return to the planning, the user must click on the

PATIENT INFORMATION

7. Planning realization

7.1. PATIENT INFORMATION

The first step is to fill in the information about the patient and the procedure.

The user must indicate the name and surname of the patient **OR** the anonymized code related to the case, if he does not wish to reveal the identity of the patient.

This information will be automatically reported in the general insert located at the top left of the application:

<u>Name and Surname</u>	<u>Anonymized ID</u>
<p>ANONYMIZED ID :</p> <input type="text"/>	<p>ANONYMIZED ID * :</p> <input type="text" value="PatientCode"/>
<p>LAST NAME * :</p> <input type="text" value="Doe"/>	<p>LAST NAME :</p> <input type="text"/>
<p>FIRST NAME * :</p> <input type="text" value="John"/>	<p>FIRST NAME :</p> <input type="text"/>

The user can also indicate the expected date of the intervention, in the format **dd/mm/yyyy**.

Optionally, the user can indicate the weight (in kg) and height (in cm) of the patient.

He must also fill in the operated side by clicking on  for a right hip or on  for a left hip.

This information will be automatically reported in the general insert located at the top left of the application:

<u>Left side</u>	<u>Right side</u>
 ANONYMIZED ID LAST NAME : FIRST NAME :	 ANONYMIZED ID LAST NAME : FIRST NAME :

The user must indicate whether to perform a planning in free mode or in advanced mode, to determine the future interactions he will have with the device:


- **Advanced mode:** The surgeon will have to draw in the "Drawing" insert a trapezoid to simulate the diaphyseal axis of the stem as well as a circle to simulate the diameter of the cup. This will automatically place the implants according to the drawings made. He


PATIENT INFORMATION

will also be able to modify the positioning and orientation of his implants with dedicated tools.

- Free mode: The surgeon will have the implants randomly placed on the layer. He will have to carry out all the positioning of the implants himself.


Finally, the surgeon must load onto the planner the radio on which he wishes to carry out his planning. To do this, it can:

- Click on . A browser window opens and the user can browse their folders to select imaging.
- Click in the drop area directly. A browser window opens and the user can browse their folders to select imaging.
- Drag and drop by first choosing imaging from a folder and dragging it into the area provided for this purpose on the planner :


In both cases, imaging will be displayed on the planner:

PATIENT INFORMATION


The user must then fill in one last piece of information on this page:

If the loaded image contains only one hip, he must come and click on

If the loaded image contains the entire pelvis, he must come and click on


Once all the information mentioned above has been filled in, the user can then click on the


which automatically appears at the bottom right of the screen or by

clicking on the navigation arrow .

A confirmation message is displayed:


By clicking on , The surgeon validates the patient's consent for the use of his data for statistical processing in order to design medical devices based on artificial intelligence models. Then, the CALIBRATION page will be active and visible.


CALIBRATION

7.2. CALIBRATION

Once the "Calibration" navigation arrow is activated  , A new page is displayed on the screen:


If the automatic calibration worked, the user will find a validation message :

He then finds the result of this automatic calibration, schematized by a circle automatically superimposed on the ball present on imaging :

CALIBRATION


Repositioning the calibration circle by clicking on the orange area

Resizing the diameter of the calibration circle

The user can change the calibration performed automatically by moving the resizing point or travel area described above. By hovering the mouse over the inside of the circle, the  will be displayed and the user can then delete the displayed circle. He then has the opportunity to redraw a new circle.

If the automatic calibration did not work, an alert message will be displayed on the screen :

The user will have access to the same tools to size and position himself the circle allowing the calibration of imaging.

On this page, the value of the calibration ball indicated by default in the surgical preferences if existing. This value can be modified in the field found in the explanatory menu to the left of

imaging :

CALIBRATION BALL (mm) :	28
--------------------------------	-----------

Under imaging, the user will find sliders for adjusting the contrast and brightness of the loaded radio .

CALIBRATION

Once the calibration is done, depending on the type of schedule chosen, the user can then choose between 2 options:

- For **Advanced Planning**,

- Click on the navigation arrow


- Click on the navigation arrow

In this case, the user can go to paragraph 7.3 of this document to continue his planning.

- For **Free Planning**,

- Click on the navigation arrow

- Click on the navigation arrow


In this case, the user can go to paragraph 7.4 of this document to continue his planning.

DRAWING

7.3. DRAWING


Once the "DRAWING" navigation arrow is activated, a new page is displayed on the screen:

By default, the **CIRCLE** button is on. The first action is therefore to draw a circle to simulate the diameter of the femoral head. The center of the head will be used for the automatic positioning of the cup in the "IMPLANTS" page. For this, the user will select with a left click of the mouse the center of the head and, without releasing the click of the mouse, he will stretch the circle by dragging his mouse:

Repositioning of the circle simulating the femoral head

Resizing the diameter of the circle simulating the femoral head


As with calibration, the user can adjust the drawing made by moving the sizing point and clicking in the green area to move the circle, as shown above. The user can also reset his drawing to start over by clicking .

DRAWING

The user must then activate the tool **TRAPEZOID** by clicking on it. He can then come and click the 4 points representing the 2 levels of cutting of the femoral diaphysis. The midpoints of these bases make it possible to automatically define the femoral diaphyseal axis:

Repositioning of the simulated diaphyseal axis

Resizing trapezoid points

The user can adjust the drawing made by moving the points described above. The user can also reset his drawing to start over by clicking  in the center of the trapezoid.

Under imaging, the user will find sliders for adjusting the contrast and brightness of the loaded radio .

Once all the drawings mentioned above have been made, the user can then click on the navigation arrow  or  to go to the next page.

IMPLANTS


7.4. IMPLANTS

Once the "IMPLANTS" navigation arrow is activated


, a new page is displayed on the screen :

For Advanced Planning :


For Free Planning :

IMPLANTS

Depending on the planning mode chosen, the user can perform different movements on the stem and on the cup:

- Stem position button on layer:
 - Advanced planning: 3 types of possible trips:
 - Rotate the stem around the Z-axis
 - Free translation in X and Y
 - Translation along the diaphyseal axis
 - Free planning:
 - Rotate the stem around the Z-axis
 - Free translation in X and Y

- Button to position the **cup** on the layer:
 - Advanced and free planning: 2 types of possible trips:
 - Rotation of the cup around the Z axis
 - Free translation in X and Y

On the layer, the tilt of the cup is displayed. The default value is 45°. This value is modified according to the orientation given by the user via the tools available.

Regardless of the planning mode chosen, the user also has access to various configuration, sizing and visualization tools :

- For the cup:

IMPLANTS

The screenshot shows the ONEORTHO software interface for selecting a cup. At the top, it displays "ONEORTHO" and "Cementless dual mobility". Below this, there is a "CUP SIZE" section with a minus sign, a green box containing the number "52", and a plus sign. To the left of the cup size is a small image of a cup. Below the cup size, there are two color selection options: "Color:" with a green circle and a white circle, and "Other:" with a green circle. Below these are two toggle switches: "Plan cup on opposite side:" and "Display cup on opposite side:". A callout box on the left explains that adding a cup on the opposite side to the intervention (if this button is active, the above tools will be active for the opposite side). A callout box on the right explains that the display of a cup on the side opposite the intervention (if this button is active and the planning button is inactive, the above tools will be active for the operated side). Another callout box on the right explains that the cup color selection (preset or free choice) is controlled by the "Other:" option.

Choice of cup type

Choice of cup size

Adding a cup on the opposite side to the intervention (if this button is active, the above tools will be active for the opposite side)

Cup color selection (preset or free choice)

Display of a cup on the side opposite the intervention (if this button is active and the planning button is inactive, the above tools will be active for the operated side)

IMPLANTS

- For the stem:

Choice of stem type

Choice of stem neck type

Choice of stem size

Stem color selection (preset or free choice)

Display of a stem on the side opposite the intervention (if this button is active and the schedule button is inactive, the above tools will be active for the Operated)


Adding a rod on the opposite side of the procedure (if this button is active, the above tools will be active for the opposite side)

The user can also observe in real time the value of the prosthetic length and offset :

An diagram appears on the screen when a tap  is done :

IMPLANTS

Under imaging, the user will find sliders for adjusting the contrast and brightness of the loaded radio .

For Advanced planning only, it can also affect the opacity of the circle and trapezoid drawings on the layer (Displayed full, hidden at 50%, hidden at 100%) :

Once the stem is planned, the user can simulate the cut to display the prosthetic head. To


do this, the user must click on  or  to go to the next page.

CUTS

7.5. CUTS


Once the "CUTS" navigation arrow is activated, a new page is displayed on the screen:

The objective of this page is to be able to carry out the simulation of the femoral cut and therefore to determine the diameter of the head of the stem as well as the length of its neck.

To do this, the user must come and click on the imagery the 2 points tracing the right of the cutting plane of the femoral stem. These two points are movable, by clicking on them and holding down to move. Releasing the mouse validates the movement. 3 other points, also editable, is displayed on the screen:


CUTS

Once the points are positioned, the user can start the simulation by clicking on

in the side menu to the left of the planning. The cut simulation is then displayed on the screen and the button becomes grayed out :


If the user wants to change the points of the cutting simulation and restart it, click on

The "Simulate cut" button will then be active again.

In the side menu to the left of the planning, the user also finds the tools specific to the sizing of the diameter of the prosthetic head and the length of the neck.:

The modification of these values causes the simulation of the cut to be modified simultaneously on imaging.


VALIDATION

7.6. VALIDATION AND REPORTING


Once all the implants are sized and positioned, the surgeon validates his planning by clicking

on . A message appears on the screen :

By clicking on , The surgeon can return to the different tabs proposed and modify his planning.

By clicking on , the surgeon validates his planning and a report gathering the important elements of it is generated and downloaded automatically (see Annex 1)

The user will be redirected to archived planning.

MAINTENANCE - SECURITY

8. Maintenance

An update is carried out annually to carry out security maintenance, it is automatic and completely transparent for the user. This maintenance includes resolving bugs and updating dependencies such as the framework, function libraries and SOUPs. It aims to prevent the emergence of new bugs and security vulnerabilities.

Depending on their criticality, some bugs can be resolved without delay.

9. Cybersecurity

Our approach to mastering cybersecurity is based on risk management at each stage of the software lifecycle in order to integrate cybersecurity concepts by design. Security measures are applied regarding user authentication to the web application which is managed by login and password that allows to give access only to a user to his own cases and therefore to restrict the access of patient data to authorized persons only. For this we ask you to be careful not to communicate your password and to log out when you stop using the application.

All risks related to cybersecurity are integrated and controlled in the risk management file of the medical device.

10. Contact

Contact :

For any questions concerning the use of the platform, please contact ONE ORTHO at the following email address: contact@oneortho-medical.com or by phone at 04 26 78 76 74.

The ONE ORTHO team is available to answer you from Monday to Friday from 9am to 6pm.

All risks are controlled through our risk analysis and no residual risk is present. Despite the rigour we put into the development of this software, any serious incident in connection with the device must be notified to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

A list of known bugs is available at the following link:

<https://github.com/oneorthomedical/2D-HIP-PLANNING>

In the event that a bug not present in this list is detected, please send the information to the ONE ORTHO service.

MAINTENANCE - SECURITY

11. Manufacturer

ONE ORTHO
 Parc INOPOLIS,
 206 Route de Vourles
 69230 Saint Genis Laval
 FRANCE
 contact@oneortho-medical.com

Class IIa medical device, manufactured by ONEORTHO Medical, complies with the requirements of EU Regulation 2017 / 745.

Date of obtaining the 1st CE marking: 2017.

© 2023 OneOrtho Medical. All rights reserved.

12. Information concerning the device and associated symbols

Light Theme :

 Manufacturer	 Date of manufacture	 Unique device identifier	 Medical device
 Consult instructions for use or consult electronic instructions for use			

Dark Theme :

 Manufacturer	 Date of manufacture	 Unique device identifier	 Medical device
 Consult instructions for use or consult electronic instructions for use			

MAINTENANCE - SECURITY

Copyright Information

This document and what it describes are protected by the Law of 11 March 1957 on Literary and Artistic Property, amended by the Law of 3 July 1985, as well as by copyright laws and international conventions.

Without the written permission of OneOrtho Medical no part of this manual, including the products and software described herein, may be reproduced, transcribed, stored in a database system, or translated into any language, form or means, except documentation retained by Buyer for backup purposes. Any person who does not comply with these provisions will be guilty of counterfeiting and will be liable to criminal penalties provided for by law.

OneOrtho Medical provides this manual "as is" without warranty of any kind, either express or implied, or on terms of merchantability or fitness for a particular purpose. In no event shall OneOrtho Medical, its directors, officers, employees, agents be liable for any indirect, special, incidental or consequential damages even if OneOrtho Medical has been advised of the possibility of such damages arising from any defect or error in this manual or product.

The names of products and companies that appear in this manual may or may not be registered trademarks, or subject to copyright for their respective companies, and are used only for identification or explanation purposes, and for the sole benefit of the owners, without intent to infringe.

ANNEX

13. ANNEX

13.1. Annex 1 : Planning report

PLANNING REPORT — 2D HIP —

Case number
Side

PATIENT INFORMATION

Anonymous code / LAST

NAME

First name

Case number

SURGERY INFORMATION

Surgery date


Surgeon

Planning date & time

2D PLANNING SUMMARY


	CUP	LINER	STEM	HEAD
DESCRIPTION / REFERENCE				
SIZE				

OFFSET		LENGTH	
--------	--	--------	--


IMPLANT RESIZING COEFFICIENT:		BALL SIZE:	
-------------------------------	--	------------	--

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