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
The screenshot displays the OneOrtho medical software interface for a preoperative knee assessment. The main view shows a 3D model of the knee joint with various measurements and callouts. The interface includes a top navigation bar with 'Preoperative View' and 'Planning' tabs, and a right sidebar with patient information: 'Knee left', 'Surgical Ref: DEMO-Left', 'Last Name: DEMO', 'First Name: Left', 'Born on: 01/01/1976', and 'Surgery: 01/06/2026'. A warning message states: 'This view is based on the acquisition of the patient's leg in the imaging system. The joint is not weight-bearing.' The measurements shown are:

- HKS 6.2°**: Angle between the femoral mechanical axis and the femoral anatomical axis in the frontal plane. Measured internally.
- LDFA 86.8°**: Angle between the femoral mechanical axis and the axis passing through the distal femoral condyles (medial and lateral) in the frontal plane. Measured externally.
- HKA 176.2°**: Angle between the femoral mechanical axis and tibial mechanical axis in the frontal plane. Measured internally.
- MPTA (Varus) 87.2°**: Angle between the tibial mechanical axis and the axis passing through the medial and lateral centers of the tibial plateau, internally, in the frontal plane.
- Flexion 10.1°**: Angle between the tibial mechanical axis and the femoral mechanical axis in the sagittal plane (complementary to the large angle). Flexion: knees in flexion. Recurvatum: knees in extension.
- Tibial slope 13.1°**: Angle between the perpendicular to the anteroposterior median axis and the tibial mechanical axis in the sagittal plane.
- Internal space 1.6mm**: Distance between the distal end of the internal part of the condyle and the bottom of the internal tibial cup.
- External space 5.5mm**: Distance between the distal end of the external part of the condyle and the bottom of the external tibial cup.
- 833 mm**: Sum of the lengths of the femoral mechanical axis and the tibial mechanical axis in 3D (non-projected).

Buttons for 'Show all info' and 'Hide all info' are visible at the bottom of the interface.

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
The screenshot displays the OneOrtho software interface for knee planning. The interface is divided into several sections:

- Top Bar:** Includes patient information (Knee left, Surgical Ref: DEMO-Left), last name (DEMO), first name (Left), and surgery date (02/19/2026).
- Left Panel:** Contains three sections for bone planning:
 - EMUR (Patella):** Reference: Anterior, Rot. axis: Distal femoral axis, Meas. ref.: Posterior condyles, Range: CONDYLE PS, Size: 7.
 - PATELLA:** Reference: Posterior, Meas. ref.: Posterior condyles, Diameter: 30mm, Thickness: 10mm.
 - TIBIA:** Reference: Posterior plateau ed..., Meas. ref.: Posterior plateau ed..., Range: EMBASE FIXE, Size: 6, Insert: 5/7.
- Central 3D View:** Shows a 3D model of the knee joint with various views (Anterior, Lateral, etc.) and measurement tools. The 'External rotation', 'Posterior resection', and 'Flexion' controls are highlighted in purple. The 'Views' button is highlighted in green.
- Right Panel:** Contains a 'Views' button (highlighted in green) and a 'Colors' button.

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
OneOrtho LEFT Genu Verso 3.8" Cartilage not included in CT imaging Initial HGA 176.2" Preoperative View Planning INFO PATIENT Knee left Surgical Ref: DEMO-Left Last Name: DEMO First Name: Left Born on: 01/01/1976 Surgery: 02/19/2026

PATELLA
 Reference: Posterior
 Meas. ref: Posterior condyles
 Diameter: 30mm
 Thickness: 10mm

External rotation: 3°
Anterior notching: 1.5 mm
Femoral varus: 0.5°
Distal resection: 8.4 mm
Flexion: 3.5°

Views (highlighted)

Planification os rotulien:

OneOrtho LEFT Genu Verso 3.8" Cartilage not included in CT imaging Initial HGA 176.2" Preoperative View Planning INFO PATIENT Knee left Surgical Ref: DEMO-Left Last Name: DEMO First Name: Left Born on: 01/01/1976 Surgery: 02/19/2026

External rotation: 0°
Posterior resection: 8 mm
Flexion: 0°

Views (highlighted)

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The screenshot displays the OneOrtho software interface for knee surgery planning. The interface is divided into several sections:

- Top Bar:** Includes the OneOrtho logo, patient information (LEFT, Genu Varum 3.8°, Cartilage not included in CT imaging, Initial HKA 175.2°), and patient details (Knee left, Surgical Ref: DEMO Left, Last Name: DEMO, First Name: Left, Born on: 01/01/1976, Surgery: 02/19/2026).
- Navigation and Tools:** A central toolbar with icons for actions like edit, copy, paste, and navigation.
- Parameter Panels:**
 - FEMUR:** Reference: Anterior; Rot. axis: Diaphyseal axis; Meas. ref.: Posterior condyles; Range: CONDYLE PS; Size: 7.
 - TIBIA:** Reference: Posterior plateau ed...; Meas. ref.: Posterior plateau ed...; Range: EMBASE FIXE; Size: 6; Insert: 6/7.
- 3D Views and Measurements:** Multiple views of the knee joint with various measurements:
 - External rotation:** 3°
 - Anterior notching:** 1.5 mm
 - Femoral varus:** 0.5°
 - Distal resection:** 8.4 mm
 - Flexion:** 3.5°
 - External rotation (bottom):** 0°
 - Genu Varum:** 0°
 - Tibial resection:** 9.5 mm
 - Tibial slope:** 3°
- Dimensions:** Various linear measurements are shown, such as 11.8 mm, 20.3 mm, 8.5 mm, 9 mm, 18.5 mm, 9.5 mm, 8.4 mm, 16.9 mm, 8.5 mm, 179.6°, 4.8 mm, 14.3 mm, 9.5 mm, 2.2°, and 0.6 mm.
- Views:** A vertical toolbar on the right includes icons for 3D, Views, and Colored cuts.

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
oneOrtho LEFT Genu Varum 3.0°
Cartilage not included in CT imaging
Initial HGA 176.2°

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OneOrtho med fcs | LEFT | Genu Varus 3.P | Cartilage not included in CT imaging Initial HKA 176.2°

Preoperative View | Planning

INFO PATIENT | Knee left | Surgical Ref: DEMO-Left | Last Name: DEMO | First Name: Left | Born on: 01/01/1976 | Surgery: 02/19/2026

External rotation: 9° | Anterior notching: 1.5 mm | Femoral varus: 0.5° | Distal resection: 8.4 mm | Flexion: 3.5°

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The image shows a software interface for planning a patella implant. At the top, there is a green button labeled 'PATELLA'. Below it, there are five rows of configuration options:

- Reference:** A dropdown menu with 'Posterior' selected.
- Meas. ref.:** A dropdown menu with 'Posterior condyles' selected.
- Range:** A dropdown menu with 'RESURFACAGE' selected.
- Diameter:** A numeric input field with '34mm' and green minus/plus buttons.
- Thickness:** A numeric input field with '8mm' and green minus/plus buttons.

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www.oneortho-medical.com

TABLE OF CONTENTS:

1-	<u>GENERAL</u>	77
2-	<u>TERMS OF USE</u>	83
3-	<u>LIFETIME OF THE DEVICE</u>	84
4-	<u>GENERAL INFORMATION ON DISPLAYS</u>	85
5-	<u>SURGICAL PREFERENCES</u>	89
6-	<u>“PREOPERATIVE VIEW” PAGE</u>	90
7-	<u>PLANNING</u>	95
	7.1. MODIFICATION OF PLANNING	95
	7.1.1. Planning the femur	111
	7.1.2. Planning the tibia	116
	7.1.3. Planning the patella	120
	7.1.4. 3D planning visualization	123
	7.1.5. Indication of space generated	125
	7.2. PLANNING VALIDATION AND REPORT GENERATION	127
8-	<u>MAINTENANCE</u>	129
9-	<u>CYBERSECURITY</u>	129
10-	<u>INCIDENT / BUG</u>	129
11-	<u>CONTACT</u>	129
12-	<u>MANUFACTURER</u>	130
13-	<u>DEVICE INFORMATION AND ASSOCIATED SYMBOLS</u>	130
14-	<u>APPENDIX 1 : GETTING STARTED GUIDE</u>	131

Restrictions on use:

This software is intended for use by professionals (Orthopaedic surgeons). In no way should it be used as the sole basis for clinical decisions on diagnosis, care or management of the patient. The plausibility of information obtained through the software must be systematically verified clinically before it is used to treat patients. Any application of the program's medical information that would not be in line with the original concept or intended use of the program is not recommended and will be considered as misuse of the software.

General

1- General

DEVICE NAME / TRADE NAME:

The name of the device range is "3D KNEE PLANNING SOFTWARE".

The reference presented in this document is the "Digital Planner Total" reference.

INTENDED USE:

3D KNEE PLANNING SOFTWARE is intended to provide orthopaedic surgeons with basic support for clinical decision-making through 3D visualization and planning tools for preoperative surgery planning in patients undergoing knee arthroplasty.

GENERAL DESCRIPTION:

The 3D KNEE PLANNING SOFTWARE is a planning software designed for surgical planning of knee prosthesis. It is a device that provides information to clinical management via planning tools and 3D visualization of proximal tibial, distal femoral and patellar anatomy. Mainly by that mean, it allows the planning of the positioning, size and model of total knee prosthesis in order to plan surgery alone or to plan the surgery and design a custom-made guide, if required, through the associated workflow platform. Thus, it provides the surgeon user with information that constitutes basic level of support to take a decision regarding the pre-operative planning of the surgery. It includes a 3D library with several implant models, in their different sizes, of total knee prostheses. Once the planning has been performed and validated by the surgeon user, the software provides a planning report

Using 3D items reconstructed from the DICOM images provided by the surgeon, this software displays the patient's knee (proximal tibia, distal femoral and patellar bone structure) with a 3D view and 2D views (NIfTI images).

It displays the knee prosthesis within the digitally reconstructed bone structure. The surgeon can virtually manipulate the prosthesis components (tibial, femoral, and patellar) with precision by adjusting parameters such as lateral position, femoral tibial and patellar cut depth, inclination angles, anterior notching range, and size.

The default position of the prosthesis is set by One Ortho's Digital Production Department in such a way as to reduce the case's processing time for the surgeon user and is not a planning suggestion.

SOFTWARE VERSION:

The active version of the "Digital Planner Total" software is **version 2.1**.

General

PRECAUTIONS FOR USE:

Residual risks	Recommendations	Patient consequences
<p>The user does not understand how to open or use the 3D scene.</p> <p>The user is unable to find a value's definition using the "?" button which appears when hovering the value.</p> <p>The user cannot find the online help section. (usability risks)</p>	<p>Before using the device, go through the User Manual: all steps and tools are described.</p> <p>Also go through the Appendix 1 "Getting started Guide".</p> <p>Always follow good practices related to the surgery, see section "Risk control measures external to the software" stated below.</p>	<p>Longer operating time.</p>
<p>The user leaves the workstation without logging off. (usability risk)</p> <p>The user shares his login information. (usability risk)</p> <p>Access to the software by someone who should not have access rights.</p>	<p>Never transfer access rights to a third party, login and password are strictly personal.</p> <p>Always log off after you have finished using the software.</p>	<p>No patient consequence.</p>
<p>If the design of the labeling format is not carried out according to regulatory requirements (logos and information)</p>	<p>Make sure to understand the label. In case of doubt, consult the detailed information in the User Manual.</p>	<p>No patient consequence</p>
<p>Wrong User Manual version on the IFU website leading to software use error</p> <p>The user is not informed that the User Manual has been revised.</p>	<p>Ensure that the version of the User Manual used is identical to the version identified on the software label.</p> <p>Always follow good practices related to the surgery, see section "Risk control measures external to the software" stated below.</p>	<p>Longer operating time.</p>

General

Residual risks	Recommendations	Patient consequences
<p>The user uses an outdated or unverified browser.</p> <p>A security flaw due to using an outdated browser.</p>	<p>Always use a browser compliant with the manufacturer's requirements and ensure that it is not obsolete in order to avoid security failures.</p> <p>Always follow good practices related to the surgery, see section "Risk control measures external to the software" stated below.</p>	<p>Longer operating time.</p>
<p>Software not working.</p> <p>The output data (STL files) leading to the production of custom-made cutting guides does not match the planning made by the surgeon user.</p>	<p>Follow minimum requirements concerning hardware (including screen definition) and IT networks characteristics as defined in the User Manual to ensure correct display of device's interface as to avoid inconvenience and security issues.</p> <p>Always follow good practices related to the surgery, see section "Risk control measures external to the software" stated below.</p>	<p>No patient consequence.</p> <p>The medical staff must use other means to perform the plan of the surgery.</p>
<p>Erroneous data/Display problem</p>	<p>To avoid issues with the planning output data, make sure the patient's data is correct.</p> <p>Always follow good practices related to the surgery, see section "Risk control measures external to the software" stated below.</p>	<p>Error in the data.</p> <p>Longer operating time.</p>
<p>Case canceled by mistake (usability risk)</p> <p>Loss/forgetting of planning by the surgeon</p>	<p>Make sure to validate the case to generate and archive the planning report. After validation, the planning report is available on the platform.</p>	<p>No patient consequence.</p> <p>The medical staff must use other means to perform the plan of the surgery.</p>

General

Residual risks	Recommendations	Patient consequences
Bugs, non-compliance and/or material safety issues not reported to ONE ORTHO	Please let us know if you experiment any bug (see §contact)	No patient consequence.
The user needs to use the user manual in paper form.	At the user's request, the User Manual can be supplied in paper form within 7 days of the request, at no additional cost.	No patient consequence.

Risk control measures external to the software:

The use of 3D KNEE PLANNING SOFTWARE does not replace the usual intraoperative control measures.

In order to avoid increasing operating time, the user is required to follow risk control measures external to the software:

- Verify that the data on the screen matches the patient data (name/anonymized code, date of birth, side to be operated, date of intervention...)
- When using traditional instrumentation, respect the operating technique established by the manufacturer of the prosthesis, particularly:
 - the use of traditional instrumentation dedicated to measure alignment and the needed prosthesis size.
 - the use of trial prosthesis before implanting the definitive metallic prosthesis.
- When using custom-made cutting guides, verify that the cutting planes are correct regarding the patient bony anatomy, as mentioned in the Patient-Specific Instrumentation's operating technique, by:
 - using the instrumentation to check that the cutting slots are well positioned,
 - using the alignment control rod to check the posterior tilt

The data presented by the software is not a planning suggestion, and the choice of planning is made by the surgeon user.

It is forbidden to use several web browsers or browser tabs to simultaneously access different pages of the application.

INSTALLATION SPECIFICATIONS:

Operating systems:

- Windows,
- Mac OS.

Browser:


- Google Chrome

The browser used must be compatible with WebGL technology:

<https://caniuse.com/?search=webgl>

General

The zoom level must be 100% for optimal use of the platform:

General

Screen definition:

The platform is fully responsive, i.e. it adapts to the size of the screen. This adaptation can, on small smartphone-type formats, modify the layout of the user screen. It is recommended to use a minimum of 1200px in 16:9 to display all the information in the application.

Internet connection:


An Internet connection is required to use the application. A minimum connection speed of 10 Mbps is recommended.


Equipment:

The client workstation hardware required for a comfortable user experience is:

- **Processor:** The processor must be equivalent to an Intel Core i5 (sixth-generation or higher).
- **Graphics card:** The graphics card must have a minimum video memory of 2048 MB.

ACCESSIBILITY OF DOCUMENTATION :

This User Manual can be accessed at any time by clicking on the  icon.

The regulatory label is available on the planner by hovering the mouse over the icon .

Lifetime of the device

2- Terms of use

Indications

3D KNEE PLANNING SOFTWARE is indicated to assist orthopaedic surgeons in the preoperative planning of knee arthroplasty procedures in patients with knee osteoarthritis.

Intended users

The intended users of the 3D KNEE PLANNING SOFTWARE are orthopedic surgeons with specialized knowledge in musculoskeletal anatomy, biomechanics, and reconstructive surgery, including knee arthroplasty procedures. The surgeon is responsible for reviewing and approving the preoperative surgical plan using the device.

There is no specific software training required to use the device for the first time.

Target population

3D KNEE PLANNING SOFTWARE is intended for use in male and female adult patients of advanced age who have reached bone maturity, diagnosed with knee osteoarthritis and scheduled for knee arthroplasty (age close to 65-75 years old).

The device is not intended for use in pediatric patients, pregnant and breast-feeding women, or other special populations.

Performance claims

3D KNEE PLANNING SOFTWARE provides clinical management with tools to planify knee arthroplasty surgery and predict the implant size of knee prosthesis components for individual patients undergoing knee arthroplasty. The exact femoral component size is predicted in 90% of cases, and the exact tibial component size in 75% of cases.

Clinical benefits

The clinical benefits of 3D KNEE PLANNING SOFTWARE are indirect. By enabling accurate preoperative planning, the device supports surgical workflow efficiency. Its use has been associated with a reduction in total procedure time by approximately 15% compared to conventional knee arthroplasty

Contraindications

The planning software must not be used:

- For joints other than the knee joint.
- For planning other prostheses ranges than those referenced in the library.
- When knee arthroplasty is contraindicated.

Adverse side effects

No adverse side effects have been recorded in association with the device.

Lifetime of the device

3- Lifetime of the device

This software will be maintained for two years.



During this period, One Ortho guarantees the performance of the software and undertakes to provide updates to enable its optimal use.

General information on displays

4- General information on displays

In this section, we'll describe all the information shown and the possible interactions in this insert, which is common to both pages of the application. As required, the user can navigate from the preoperative page to the planning page and vice-versa.


Whichever page is displayed, the application's top banner looks like this:

To move from one page to another, simply click on the corresponding navigation arrows:


-  The "Preoperative view" page is active.
-  The "Planning" page is active.

In the top right-hand corner of the application, the following information is displayed:


General information on displays

Below the navigation arrows, the following information is displayed:

Tools for specific actions are also available (can only be used on the Planning page):


Pressing the Colourimetry button  will display two sliders that can be used to adjust the contrast and brightness:


Pressing the button  will take a screenshot of part of the interface, which can be saved to the user's computer or added to the planning report. To do this, select a rectangular area by holding down the left mouse button and moving the mouse, while remaining within the browser window.

General information on displays


Once the button is released, the following interface is displayed to name the screenshot. Please note that if the selected area extends beyond the browser, the screenshot will be blank.

The buttons at the bottom allow you to cancel, save the image to the user's workstation, or add it to the planning report.

Pressing the button  will display the following interface, allowing you to choose the colours of the bones and implants.


Choice of bone colour

Choice of implant colour from a range of predefined colours


Free choice of implant colour

General information on displays

Tools for managing navigation in the application are available (only the temporary save button is active on the "Preoperative view" page):

Finally, tools for managing displays, accessing explanatory documentation and disconnecting from the platform are available at any time during the planning process:


Surgical preferences

5- Surgical preferences

Users can specify their preferred planning parameters.

To do this, access the preference settings screen via the  button:


The interface for this page is as follows:

The user can then set the following parameters:

- The appearance of the planning theme:

Light theme


Dark theme

Surgical preferences

- The default page displayed when the case is opened

“Preoperative view” page

“Planning” page

- **Default values for the display of cut bones (femur, tibia, and patella):**
 - Hidden
 - Displayed
 - If displayed:
 - Transparency level (0% = transparent, 100% = opaque)
 - Option to display cut surfaces in color
- **Default values for the display of uncut bones (femur, tibia, and patella):**
 - Hidden
 - Displayed
 - If displayed:
 - Transparency level (0% = transparent, 100% = opaque)

Default values for certain femoral planning information:

- Reference
- Rotation axis
- Femoral implant range

Default values for the display of the femoral implant:

- Hidden
- Displayed
- If displayed:
 - Transparency level (0% = transparent, 100% = opaque)


Default values for certain tibial planning information:

- Reference
- Tibial implant range

Default values for the display of the tibial implant:

- Hidden
- Displayed
- If displayed:
 - Transparency level (0% = transparent, 100% = opaque)


To validate any changes made and return to the planning, the user must click on the button

Preoperative view page


6- "Preoperative view" page

One of the application's pages is entitled "Preoperative view". It contains all data relating to the patient's native anatomy:

The views shown are illustrations based on the DICOM images of the patient's leg according to the established imaging protocol. The joint is not under load.


Preoperative view page

Definition of displayed measurements:

- **HKS:** Angle between the femoral mechanical axis and the femoral anatomical axis in the frontal plane. Measured internally.
- **HKA:** Angle between the femoral mechanical axis and tibial mechanical axis in the frontal plane. Measured internally.
- **LDFA:** Angle between the femoral mechanical axis and the axis passing through the distal femoral condyles (medial and lateral) in the frontal plane. Measured externally.
- **MPTA (Tibial varus):** Angle between the tibial mechanical axis and the axis passing through the medial and lateral centers of the tibial plateau, internally, in the frontal plane.
- **INTERNAL SPACE:** Distance between the distal end of the internal part of the condyle and the bottom of the internal tibial cup.
- **EXTERNAL SPACE:** Distance between the distal end of the external part of the condyle and the bottom of the external tibial cup.

Preoperative view page



Definition of displayed measurements:

- **Flexion:** Angle between the tibial mechanical axis and the femoral mechanical axis in the sagittal plane. If the value is positive, the knee is in flexion. If it is negative, the knee is in extension.
- **Tibial slope:** Angle between the perpendicular to the anteroposterior median axis and the tibial mechanical axis in the sagittal plane.

Preoperative view page

At any time, the user can view the measurement definition by clicking on it.

They also may view all the definitions by pressing the  button:

Planning

7- Planning


7.1. Modification of planning

One of the pages in the application is entitled "Planning". It contains all the tools needed to plan the current case:

Dynamic 3D visualization and 2D visualization layout insert

Interactive insert for 2D visualization, planning, movement and positioning of the femoral implant

Visualization tools access panel

Dynamic insert to indicate generated spaces

Dynamic HKA angle insert calculated in real time

Interactive insert for 2D visualization, planning, movement and positioning of the tibial implant

For patella planning, you must first activate the patella visualization (by double clicking)

The screen then becomes this one:

Planning

Interactive visualization insert for patellar planning

Interactive patellar implant planning insert

Button to press by double clicking to access to the patellar tools

The screenshot displays the OneOrtho software interface for knee implant planning. The interface is divided into several sections:

- Top Bar:** Includes patient information (Last Name: DEMO, First Name: Left, Born on: 01/01/1976, Surgery: 02/19/2026) and navigation tabs for "Preoperative View" and "Planning".
- Left Panel:** Contains three main sections: "EMUR" (knee), "PATELLA" (patella), and "TIBIA" (tibia). The "PATELLA" section is highlighted with a yellow box and contains adjustable parameters: Reference (Posterior), Meas. ref. (Posterior condyles), Diameter (30mm), and Thickness (10mm).
- Central View:** Shows a 3D model of the knee joint with various views (M, L, LATERAL) and adjustable parameters for "External rotation" (0°), "Posterior resection" (8 mm), and "Flexion" (0°). A purple box highlights these three parameters.
- Right Panel:** Contains a "Views" section with a button highlighted by a green box, and a "Colored cuts" section with a toggle switch.
- Bottom Panel:** Shows additional adjustable parameters for "External rotation" (0°), "Genu Varum" (0°), "Tibial resection" (9.5 mm), and "Tibial slope" (3°).

Planning

Description of the visualization tools access panel:

Scanner view/3D view


3D view active

Scanner view active

“3D model visualization” mode activated:

Planning

“Scanner image visualization” mode activated:

In this scanner view, the user can move section by section using the mouse wheel and view the overlaying of the selected sections in 3D mode and the selected implants. In this way, the user can see the coverage of the selected implant and validate or not his choice of planning.

A zoom slider is available for each slice in the CT scan view. The user can drag the slider upward to zoom in or downward to zoom out. The zoom is applied in real time to the corresponding slice, with no view reload required.

This feature allows the user to inspect anatomical details more closely, particularly when assessing bone structures during tibial stem extension planning.

When a slice is zoomed in, the area of interest may fall outside the visible frame. The user can navigate within the enlarged image without changing the zoom level by holding down the right mouse button on the slice and moving the mouse in the desired direction. The slice follows the mouse movement. The user releases the right mouse button to stop panning.

The zoom level remains unchanged throughout the panning operation.

Planning

Display of the Patella and planning

**Display Patella
vizzualization
tools**

**Planning of the
patella**

Patella hidden

Display Patella vizzualization tools:

The screenshot shows the OneOrtho software interface for knee implant planning. The top navigation bar includes 'Preoperative View' and 'Planning'. The patient information at the top right indicates 'Knee left', 'Surgical Ref: DEMO-Left', 'Last Name: DEMO', 'First Name: Left', and 'Born on: 01/01/1976', 'Surgery: 02/19/2026'. The main workspace displays various 3D views of the knee joint with measurement tools. The left sidebar contains settings for 'FEMUR', 'PATELLA', and 'TIBIA'. The 'PATELLA' settings panel is highlighted with a red box, showing parameters like Reference (Posterior), Meas. ref. (Posterior condyles), Diameter (30mm), and Thickness (10mm). The 'Views' toolbar on the right has the 'Patella' icon highlighted with a yellow box. The main workspace shows measurements for External rotation (3°), Anterior notching (1.5 mm), Femoral varus (0.5°), Distal resection (8.4 mm), Flexion (3.5°), and Tibial slope (3°).




Planning of the patella:

This screenshot is similar to the one above, showing the OneOrtho software interface. In this view, the 'External rotation' (0°) and 'Posterior resection' (8 mm) settings are highlighted with a red box. The 'Patella' icon in the 'Views' toolbar is highlighted with a yellow box. The 'PATELLA' settings panel is also visible, showing the same parameters as in the previous screenshot.

Patella hidden:

Planning

Visualization "Bone cut models"

Bone cut models displayed :

Bone cut models displayed transparently:

Bone cut models hidden:

Planning

OneOrtho medical LEFT Genu Varum 3.0° Certificate not included in CT imaging Initial HKA 176.2°

Preoperative View Planning

INFO PATIENT Knee left Surgical Ref: DEMO-Left Last Name: DEMO First Name: Left Born on: 01/01/1976 Surgery: 02/19/2026

Actions

ENTRANCE

FEMUR

Reference: Anterior
 Rot. axis: Diaphyseal axis
 Meas. ref.: Posterior condyles
 Range: CONDYLE PS
 Size: 7

PATELLA

Reference: Posterior
 Meas. ref.: Posterior condyles
 Diameter: 30mm
 Thickness: 10mm

TIBIA

Reference: Posterior plateau ed...
 Meas. ref.: Posterior plateau ed...
 Range: EMBASE FIXE
 Size: 6
 Insert: 6/7

External rotation: 3°
 Anterior notching: 1.5 mm
 Femoral varus: 0.5°
 Distal resection: 8.4 mm
 Flexion: 3.5°



98.0°
 0.3 mm
 20.4 mm
 8.0 mm
 11.2 mm
 6.2°
 20.4 mm
 8.0 mm
 0.0mm
 0.2° / 4.9°
 11.8 mm
 9 mm
 8.4 mm
 4.8 mm
 20.3 mm
 18.5 mm
 16.9 mm
 179.6°
 14.3 mm
 8.5 mm
 9.5 mm
 8.5 mm
 9.5 mm
 LATERAL
 2.2°

External rotation: 0°
 Genu Varum: 0°
 Tibial resection: 9.5 mm
 Tibial slope: 3°

3D
 Views
 Colored cuts

Planning

Visualization "Femoral implant"

Femoral implant displayed:


The screenshot shows the OneOrtho software interface for knee planning. The femoral implant is visible in the 3D model. The interface includes a top navigation bar with 'Preoperative View' and 'Planning' tabs. On the left, there are three panels for 'FEMUR', 'PATELLA', and 'TIBIA' with various settings. The main workspace displays multiple views of the knee joint with various measurements and adjustment sliders. A yellow box highlights the 'Femoral implant' icon in the right-hand 'Views' panel.

Femoral implant displayed transparently:

This screenshot is identical to the previous one, but the femoral implant is rendered as a transparent yellow shell. A yellow box in the 'Views' panel highlights the 'Femoral implant' icon, which now has a transparency slider next to it.

Femoral implant hidden:


Planning

Vizualization « Patellar implant »


Patellar implant displayed :


Patellar implant displayed transparently :

Planning

Patellar implant hidden :

Visualization "Tibial implant"

Tibial implant displayed

Tibial implant displayed transparently

Tibial implant hidden

Tibial implant displayed:

Planning

OneOrtho medical LEFT Genu Varum 3.8° Cartilage not included in CT imaging Initial HKA 176.2°

Preoperative View Planning

INFO PATIENT Knee left Surgical Ref: DEMO-Left Last Name: DEMO First Name: Left Born on: 01/01/1976 Surgery: 02/19/2026

Actions

3D Views

Colored cuts

Reference: Anterior
Rot. axis: Diaphyseal axis
Meas. ref.: Posterior condyles
Range: CONDYLE PS
Size: 7

Reference: Posterior
Meas. ref.: Posterior condyles
Diameter: 30mm
Thickness: 10mm

Reference: Posterior plateau ed...
Meas. ref.: Posterior plateau ed...
Range: EMBASE FIXE
Size: 6
Insert: 6/7

External rotation: 3°
Anterior notching: 1.5 mm
Femoral varus: 0.5°
Distal resection: 8.4 mm
Flexion: 3.5°

98.0°
0.3 mm 20.4 mm
8.0 mm
11.2 mm
M L
0.2° / 4.9°
11.8 mm 9 mm 8.4 mm 4.8 mm
20.3 mm 18.5 mm 16.9 mm 179.6° 14.3 mm
8.5 mm 9.5 mm 8.5 mm 9.5 mm
LATERAL
2.2°
External rotation: 0°
Genu Varum: 0°
Tibial resection: 9.5 mm
Tibial slope: 3°

Tibial implant displayed transparently:

OneOrtho medical LEFT Genu Varum 3.8° Cartilage not included in CT imaging Initial HKA 176.2°

Preoperative View Planning

INFO PATIENT Knee left Surgical Ref: DEMO-Left Last Name: DEMO First Name: Left Born on: 01/01/1976 Surgery: 02/19/2026

Actions

3D Views

Colored cuts

Reference: Anterior
Rot. axis: Diaphyseal axis
Meas. ref.: Posterior condyles
Range: CONDYLE PS
Size: 7

Reference: Posterior
Meas. ref.: Posterior condyles
Diameter: 30mm
Thickness: 10mm


Reference: Posterior plateau ed...
Meas. ref.: Posterior plateau ed...
Range: EMBASE FIXE
Size: 6
Insert: 6/7

External rotation: 3°
Anterior notching: 1.5 mm
Femoral varus: 0.5°
Distal resection: 8.4 mm
Flexion: 3.5°

98.0°
0.3 mm 20.4 mm
8.0 mm
11.2 mm
M L
0.2° / 4.9°
11.8 mm 9 mm 8.4 mm 4.8 mm
20.3 mm 18.5 mm 16.9 mm 179.6° 14.3 mm
8.5 mm 9.5 mm 8.5 mm 9.5 mm
LATERAL
2.2°
External rotation: 0°
Genu Varum: 0°
Tibial resection: 9.5 mm
Tibial slope: 3°

Planning

Tibial implant hidden :

Visualization "Bone resections"

Bone resections displayed


Bone resections displayed transparently

Bone resections hidden

Bone resections displayed:

Bone resections displayed transparently:

Planning

OneOrthro medical LEFT Genu Varum 3.8" Cartilage not included in CT imaging total HKA 176.2°

Preoperative View Planning

INFO PATIENT Knee left Surgical Ref: DEMO-Left Last Name: DEMO First Name: Left Born on: 01/01/1976 Surgery: 02/19/2026

Actions

FEMUR

Reference: Anterior
 Rot. axis: Diaphyseal axis
 Meas. ref.: Posterior condyles
 Range: CONDYLE PS
 Size: 7

PATELLA

Reference: Posterior
 Meas. ref.: Posterior condyles
 Diameter: 30mm
 Thickness: 10mm

TIBIA

Reference: Posterior plateau ed...
 Meas. ref.: Posterior plateau ed...
 Range: EMBASE FIXE
 Size: 6
 Insert: 5/7

External rotation: 3°
 Anterior notching: 1.5 mm
 Femoral varus: 0.5°
 Distal resection: 8.4 mm
 Flexion: 3.5°

98.0° 0.3 mm 20.4 mm 6.2° 20.4 mm 0.0 mm
 11.2 mm 8.0 mm 11.2 mm
 M L LATERAL
 0.2° / 4.9° 11.8 mm 9 mm 8.4 mm 4.8 mm 8.4°
 20.3 mm 18.5 mm 16.9 mm 179.6° 14.3 mm
 8.5 mm 9.5 mm 8.5 mm 9.5 mm
 2.2°

External rotation: 0°
 Genu Varum: 0°
 Tibial resection: 9.5 mm
 Tibial slope: 3°

3D Views

Colored cuts

Bone resections hidden:

OneOrthro medical LEFT Genu Varum 3.8" Cartilage not included in CT imaging total HKA 176.2°

Preoperative View Planning

INFO PATIENT Knee left Surgical Ref: DEMO-Left Last Name: DEMO First Name: Left Born on: 01/01/1976 Surgery: 02/19/2026

Actions

FEMUR

Reference: Anterior
 Rot. axis: Diaphyseal axis
 Meas. ref.: Posterior condyles
 Range: CONDYLE PS
 Size: 7

PATELLA

Reference: Posterior
 Meas. ref.: Posterior condyles
 Diameter: 30mm
 Thickness: 10mm

TIBIA

Reference: Posterior plateau ed...
 Meas. ref.: Posterior plateau ed...
 Range: EMBASE FIXE
 Size: 6
 Insert: 5/7

External rotation: 3°
 Anterior notching: 1.5 mm
 Femoral varus: 0.5°
 Distal resection: 8.4 mm
 Flexion: 3.5°

98.0° 0.3 mm 20.4 mm 6.2° 20.4 mm 0.0 mm
 11.2 mm 8.0 mm 11.2 mm
 M L LATERAL
 0.2° / 4.9° 11.8 mm 9 mm 8.4 mm 4.8 mm 8.4°
 20.3 mm 18.5 mm 16.9 mm 179.6° 14.3 mm
 8.5 mm 9.5 mm 8.5 mm 9.5 mm
 2.2°

External rotation: 0°
 Genu Varum: 0°
 Tibial resection: 9.5 mm
 Tibial slope: 3°

3D Views

Colored cuts

Planning

Change in sagittal view orientation


Medial sagittal view (left case):

The screenshot shows the OneOrtho software interface for a knee joint planning session. The main view displays a 3D model of the knee joint in a medial sagittal orientation. The interface includes a top navigation bar with 'Preoperative View' and 'Planning' tabs. On the left, there are three panels for 'FEMUR', 'PATELLA', and 'TIBIA', each with various parameters like reference, range, and size. The central area contains several adjustment sliders for 'External rotation', 'Anterior notching', 'Femoral varus', 'Distal resection', and 'Flexion'. Below these are more sliders for 'External rotation', 'Genu Varum', 'Tibial resection', and 'Tibial slope'. The right side features a vertical toolbar with various icons, including a green-bordered button with a double-headed arrow. The bottom right corner has a 'Colored cuts' toggle.

Lateral sagittal view (left case):

The screenshot shows the OneOrtho software interface for a knee joint planning session, similar to the previous one but in a lateral sagittal orientation. The main view displays a 3D model of the knee joint in a lateral sagittal orientation. The interface elements are identical to the previous screenshot, including the top navigation bar, left-side parameter panels, central adjustment sliders, and right-side toolbar. The green-bordered button with a double-headed arrow is highlighted in the toolbar.

Displaying/hiding axes


Axes displayed


Axes hidden

Planning


Axes displayed:


Axes hidden:

Enabling/Disabling coloured cuts


Colourful cuts disabled


Colourful cuts displayed

Colourful cuts displayed :

Planning


Colourful cuts disabled :

7.1.1. Planning the femur

The user is provided with a range of tools for planning the femur. They are identified on the planner by the "Purple" color. The tools for positioning and sizing the femoral implant are as follows:

Planning


The user can then select the following parameters:

<u>Reference</u>	The choice of this reference will define the default anteroposterior (AP) positioning of the condyle.
<u>Rotation axis</u>	The choice of the axis of external rotation will determine how the condyle rotates in the horizontal plane.
<u>Measurement reference</u>	The choice of this reference will define the 0 value for external rotation according to the chosen reference frame.
<u>Femoral implant range</u>	One of the distributor's ranges
<u>Size</u>	Guaranteed size in the ranges offered by the distributor


Planning

The following tools are available for making femoral cuts:

The user can then modify all these parameters to plan the femoral cuts for their case.

The indicative measurements displayed induced by the cut settings are as follows:


Planning


The display of femoral views can be modified according to the user's wishes. To do so, click on

to the right of the planner:

This will open a hidden window containing the 3D planning view and the settings for the various views:

Visualization of the distal femoral view

Planning


Visualization of the frontal femoral view

Visualization of the sagittal femoral view


These views can be combined by clicking one after the other.

For example, the user can choose to view only the distal and frontal femoral views:


The menus for sizing/positioning the femoral implant and creating femoral sections are displayed with the views selected by the user. This view is closed by clicking again on the button.

In distal femoral view only, the user can position the implant on the bone model. Passing the mouse over the femoral implant in distal view will change its appearance:

By clicking on the implant in this view, the user can move it medio-laterally. The implant then takes on a new appearance during this move:

Planning

This displacement is visible in all femoral views.

7.1.2. Planning the tibia

The user is provided with a range of tools for planning the tibia.

They are identified on the planner by the "Blue" color.

Tools for positioning and sizing the tibial implant include:

TIBIA


Reference	Posterior plateau ed... ▾
Meas. ref.	Posterior plateau ed... ▾
Range	EMBASE FIXE ▾
Size	<div style="display: flex; align-items: center; gap: 5px;"> - 6 + </div>
Insert	6/7

The user can then select the following parameters:

<u>Reference</u>	The choice of this reference will define the default positioning of the tibial component in rotation on the tibial section.
<u>Measurement reference</u>	The choice of this reference will define the 0 value of the tibial component rotation according to the chosen reference.
<u>Tibial implant range</u>	One of the distributor's ranges
<u>Size</u>	Guaranteed size within the range offered by the distributor (+/- one size of that chosen for the condyle)
<u>Insert</u>	This is an indication of the plastic spacer size compatible with the selected size.

Planning

The tools for making tibial cuts are as follows:


The user can then modify all these parameters to plan the tibial cuts for their case.

The indicative measurements displayed induced by the cut settings are as follows:

Proximal tibial sections


Planning

Proximal tibial angle


The display of tibial views can be modified according to the user's wishes. To do so, to the

right of the planner, click on >> :

This will open a hidden window containing the 3D planning view and the settings for the various views:

Planning

Visualization of the proximal tibial view


Visualization of the frontal tibial view

Visualization of the sagittal tibial view

These views can be combined by clicking one after the other.
For example, the user can choose to view only the proximal and frontal tibial views:

Planning


The menus for tibial implant sizing/positioning and tibial sectioning are displayed according to the views selected by the user. This set view is closed by clicking again on the button .

In the distal tibial view only, the user can position the implant on the bone model. Passing the mouse over the tibial implant in distal view will change its appearance:

By clicking on the implant in this view, the user can move it medio-laterally. The implant then takes on a new appearance during this move:

This displacement is visible in all tibial views.

7.1.3. Planning the patella

The user has in his possession various tools in order to carry out the planning of the patella.

Planning

They are identified on the planner by the color "Green".

The tools related to the positioning and sizing of the patellar implant are as follows:

The user will then be able to choose the following parameters:


<u>Reference</u>	The choice of this reference will define the default anterior/posterior positioning of the patellar implant.
<u>Measurement Reference</u>	The choice of this reference will set the value 0 of the external rotation according to the chosen reference. This reference is common to that of the femur.
<u>Patellar implant range</u>	One of the distributor's ranges
<u>Diameter</u>	The choice of the diameter of the patellar implant from the available sizes offered by the distributor.
<u>Thickness</u>	The choice of patellar implant thickness from the available sizes offered by the distributor.

Planning

The tools for making patellar cuts are as follows:

Rotation of the patellar implant
in the axial plane

Rotation of the patellar implant
in the sagittal plane


The user will then be able to modify all these parameters to plan the patellar resections of his case.

The indicative measurements displayed induced by the adjustment of the cuts are as follows:

Patellar offset : distance between the
prosthetic patellar apex and the
prosthetic trochlear floor

Thickness of the patella after
resection and thickness of the
resection


Remaining patellar thickness: in extension, the
distance between patellar resection and
femoral resection

Planning


7.1.4. 3D planning visualization

Throughout the planning process, the user can view the implant sizing and position on the 3D bone models of the current case. To do so, click on the  button on the right-hand side of the planner:

Planning


The following 3D view is displayed:

Users can move the entire 3D model by holding down the left mouse button on the scene presenting the femur, tibia and patella. They can also use the following visualization tools:

Visualization of bone models in frontal view

Visualization of bone models in sagittal view

Planning

Display/Hide the femur

Display/Hide the tibia

Display/Hide the patella


7.1.5. Indication of space generated

On the planning page, the user has access to the space values generated by the cuts:


The sum of the femoral and tibial sections gives the space generated by the sections, and thus the implant thickness to be added. Each measurement is the sum of the internal or external cuts, in two different configurations: flexion or extension. The four sums obtained are as follows (illustrated below):

- **Internal space in flexion:** Internal posterior femoral section + Internal proximal tibial section
- **External space in flexion:** External posterior femoral section+ External proximal tibial section
- **Internal space in extension:** distal medial femoral section + proximal medial tibial section
- **External space in extension:** distal external femoral section+ proximal external tibial section

Planning

Planning


7.2. Planning validation and report generation

Once the surgeon has set up all their cuts, dimensioned and positioned the implants, the user must validate definitively his planning. As a result, the user will no longer be able to modify what he been done, and the validated output data will be used by ONE ORTHO to produce the associated custom guides.

To do this, the surgeon clicks on the  button in the upper menu of the planning screen. They must then confirm the validation of their case by clicking on "Yes" in the window that appears:

Planning

Once validated, a report summarizing the planning carried out will be issued to the surgeon. The user is redirected to the procedure follow-up.

Maintenance, Cybersecurity & Contact

8- Maintenance

A yearly update is carried out for safety maintenance, automatically and completely transparently to the user.

9- Cybersecurity

Our approach to managing cybersecurity is based on risk management at every stage of the software lifecycle, integrating cybersecurity concepts right from the design stage. Security measures are applied to user authentication to the web application, which is managed by a login and password that gives only one user access to his or her own cases, thus restricting access to patient data to authorized persons only. We therefore ask you to be careful not to divulge your password, and to log out when you stop using the application.

All cybersecurity-related risks are integrated and controlled in the medical device's risk management file.

10- Incident / bug

Despite the rigor with which we have developed this software, any serious incident involving the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

A list of known bugs is available at the following link :

<https://github.com/oneorthomedical/3D-KNEE-PLANNING>

Should a bug not listed here be detected, please forward the information to ONE ORTHO.

11- Contact

Contact :

If you have any questions about how to use the platform, please contact ONE ORTHO at contact@oneortho-medical.com or +33 4 26 78 76 74.

The ONE ORTHO team is available to answer your questions from Monday to Friday, 9am to 6pm.

Maintenance, Cybersecurity & Contact

12- Manufacturer

ONE ORTHO
Parc INOPOLIS,
206 Route de Vourles
69230 SAINT-GENIS-LAVAL
FRANCE
contact@oneortho-medical.com

OneOrtho
medical

Class IIa medical device, manufactured by ONE ORTHO, compliant with the requirements of EU Regulation 2017/745.

Year of 1st CE marking (under EU Regulation 2017/745): 2026
2025 One Ortho. All rights reserved.

13- Device information and associated symbols

Light theme:

 Manufacturer	 Date of manufacture	 Unique device identifier	 Medical device
 Consult the User Manual			

Dark theme:

 Manufacturer	 Date of manufacture	 Unique device identifier	 Medical device
 Consult the User Manual			

Appendix

14- Appendix 1 : Getting started guide

ONE ORTHO MEDICAL

3D Knee Planner

Quick start guide
SURGEON PROFILE

OneOrto TutoStart

OneOrto medical

The image shows a laptop displaying the 3D Knee Planner software interface. The interface features several 3D models of knee joints in various views (anterior, posterior, lateral, medial, and distal). The models are rendered in a light green color. The software interface includes various control panels and data fields around the models. The OneOrto Medical logo is visible in the top right corner of the slide, and the OneOrto TutoStart logo is in the bottom left corner. The OneOrto medical logo is in the bottom right corner.

Content

- 01 General
- 02 Getting started: connection and interface
- 03 Essential functions

2

The image shows a large green vertical bar on the left side of the slide containing the ONE ORTHO MEDICAL logo. The logo consists of a stylized green 'U' shape with a white 'O' inside, surrounded by the text 'ONE ORTHO MEDICAL' in a circular arrangement. The text 'ONE ORTHO' is at the top and 'MEDICAL' is at the bottom. The background of the slide is a light beige color with a green curved shape in the bottom right corner. The content area is on the right side of the slide, featuring the title 'Content' and a list of three items: '01 General', '02 Getting started: connection and interface', and '03 Essential functions'. The number '2' is located in the bottom right corner of the slide.

Appendix

01. General

The 3D Knee Planner is a module designed for preoperative 3D planning of knee surgery.

Key features include:

1. Management of surgical preferences
2. Visualisation of preoperative imaging
3. Femur planning
4. Tibia planning
5. 3D visualisation of planning
6. Indications on the spaces generated
7. Planning validation and report generation

Installation requirements :

- **Operating systems:** Windows, Mac OS
- **Browser:** Google Chrome
The browser used must be compatible with WebGL technology (more info: <https://caniuse.com/?search=webgl>)
- **Screen resolution:** the platform is fully responsive, but we recommend using a minimum of 1200px in 16:9 to display all the information in the application, and setting the zoom level to 100%.
- **Internet connection:** an internet connection is required to use the application. A minimum connection speed of 10 Mbps is recommended.
- **Residual bugs:** if you encounter a problem with the planner, please consult the following link before contacting your manufacturer: <https://oneortho.atlassian.net/wiki/external/NmZmMjAxNDNkNjhmNGFhOGIwOWJkN2UyYjYzMDAyNmM>. An alternative solution will be proposed.

OneOrtho

3

02. First steps: connection & interface

OneOrtho

4

Appendix


Getting started: login & interface

TutoStart OneOrto

Login

CONNECTING TO / DISCONNECTING FROM THE PLATFORM

- 1 On the login page, enter:
 - Username: your email address
 - Password: set when you first logged in
 - Click on "Connection".
- 2 To log out, from the dashboard:
 - Click on your username
 - Click on "Logout"

OneOrto


5

Getting started: connection & interface

OneOrto TutoStart

Interface

- 1 ACTIONS TO PERFORM
Quick access to priority actions.
- 2 INDICATORS
Estimated average time taken to process a case and the main stages in the process.
- 3 BUSINESS MODULES
Quick access to the table of cases in progress, by joint.
For 3D knee planning, click on "Digital Planner".
Then click on the **PLANIFICATION** button.
- 4 NEWS
News from your manufacturer or OneOrto.
- 5 SETTINGS & DOCUMENTATION
 - Profile information
 - Security centre (Multi Factor Authentication)
 - File a claim
 - FAQ
 - Documentation to download

OneOrto

6

Appendix

03. Key features

OneOrto

7

Banner / Navigation menu

TutoStart OneOrto

Navigation arrows to move from the "Preoperative view" page to the "Planning" page. The green colour indicates the active page.

Information on the surgery

Patient information

Date of surgery

Native anatomical values specific to the case

Add a comment

Cancel previous action

Save temporarily

Validate the schedule definitively

Reminder of the lateralization of the case:
- L for a left knee
- R for a right knee

Information on how to measure the HKA angle

Take a screenshot (with the option of adding it to the report).

Reset planning (return to pre-planning values defined by the manufacturer).

Help (User manual)

Exit (Return to the tracking table)

Full screen mode

Surgical preferences

Interactive tutorial

OneOrto


8

Appendix

Surgical preferences

OneOrto TutoStart

1. Click on the  navigation bar to access the Surgical Preferences.
2. Save your preferred settings for each of your plans:
 - Light or dark display theme
 - Default page displayed when a case is opened
 - Default values for certain femoral and tibial planning information


OneOrto

9

Pre-operative views

OneOrto TutoStart

Data relating to the patient's native anatomy based on the acquisition of his leg according to the established imaging protocol. (The joint is not loaded).

Click on the measurement to view its definition.

You can find all the definitions in the User Manual by clicking on  in the navigation bar.

OneOrto

10

Appendix

Creating a plan

OneOrto TutoStart

Overall view

Dynamic insert for 3D visualisation and layout of 2D visualisations

Interactive insert for 2D visualisation, planning, movement and positioning of the tibial implant

Interactive insert for 2D visualisation, planning, movement and positioning of the femoral implant

Menu bar for accessing viewing tools

Dynamic insert indicating the spaces generated



Dynamic insert showing the HKA angle calculated in real time

OneOrto

11

Viewing tools

TutoStart OneOrto

Scanner view (DICOM) / 3D view
Possibility of moving section by section using the mouse wheel

Visualisation of the patellar bone.
3 display modes: bone displayed / planning / hidden

View the "cut bone model".
3 display modes: bone displayed / transparent / hidden


Display / hide landmarks

View the femoral implant.
3 display modes: bone displayed / transparent / hidden

Visualisation of the tibial implant.
3 display modes: bone displayed / transparent / hidden

Visualisation of bone resection (initial bone).
3 display modes: bone displayed / transparent / hidden

Change sagittal view orientation (lateral or medial)

OneOrto

12

Appendix


Planning the femur

OneOrto TutoStart

Choice of parameters

Reference	The choice of this reference will define the default anteroposterior (AP) positioning of the condyle.
Axis of rotation	The choice of the external rotation axis will define how the condyle will be rotated in the horizontal plane.
Measurement reference	The choice of this reference will define the 0 value of the external rotation according to the chosen reference.
Femoral implant range	One of the ranges offered by the distributor.
Size	Size guaranteed within the ranges proposed by the distributor.

Planning tools & Measurements

You can find full descriptions of the planning tools and measurements displayed in the User Manual on page 34 by clicking on :

in the navigation bar.

OneOrto

13

Planning the tibia

OneOrto TutoStart

Choice of parameters

Reference	The choice of this reference will define the default position of the rotational baseplate on the tibial section.
Measurement reference	Selecting this reference will define the 0 value for the rotation of the baseplate according to the chosen reference.
Tibial implant range	One of the ranges offered by the distributor.
Size	Guaranteed size within the ranges proposed by the distributor (+/- one size of that chosen for the condyle).
Insert	This is an indication of the insert size compatible with the size chosen.

Planning tools & Measurements

You can find full descriptions of the planning tools and measurements displayed in the User Manual by clicking on :

in the navigation bar.


OneOrto

14

Appendix

Planning the patella

OneOrto TutoStart

Choice of parameters

Reference	The choice of this reference will define the default anterior/posterior positioning of the patellar implant.
Measurement Reference	The choice of this reference will set the value 0 of the external rotation according to the chosen reference. This reference is common to that of the femur.
Diameter	The choice of the diameter of the patellar implant from the available sizes offered by the distributor.
Thickness	The choice of patellar implant thickness from the available sizes offered by the distributor.

Planning Tools & Measurements
 You can find all the descriptions of the planning tools and the measures displayed in the User Manual by clicking on:  in the navigation bar.


OneOrto

15

3D visualisation

& 2D view settings

OneOrto TutoStart

In the distal femoral or tibial view only, you can position/move the implant on the bone model.
 In the distal view, click on the implant to move it (it will turn red).

You can display 1 to 3 views for each zone (tibial and femoral).

Click here to open the 3D view window

OneOrto

16

Appendix


Space values generated

OneOrto TutoStart

- The sum of the femoral and tibial sections is used to obtain the space generated by the sections and thus to deduce the implant thickness to be added.
- Each measurement is the sum of the internal or external sections, in two different configurations: flexion or extension.

i More information in the User Manual by clicking on:

? in the navigation bar.

OneOrto



17

Planning validation

& report generation

OneOrto TutoStart

- Validation of the schedule is **definitive** (no changes can be made afterwards).
i The planning output will be used by OneOrto to design the associated custom guides.
- Once the plan has been validated, a **report** with all the data and values essential for the surgery is automatically generated for the surgeon.

OneOrto

18

OneOrtho
medical

For further information, please contact your
OneOrtho solution distributor.

OneOrtho - Parc Inopolis
206 route de Vouffes • 69230 Saint-Genis-Laval
Ed04-25

IMPLEMENT INNOVATIVE SOLUTIONS

OneOrtho
medical

For further information, please contact your
OneOrtho solution distributor.

OneOrtho - Parc Inopolis
206 route de Vouffes • 69230 Saint-Genis-Laval
Ed04-25

IMPLEMENT INNOVATIVE SOLUTIONS

OneOrtho
medical

For further information, please contact your
OneOrtho solution distributor.

OneOrtho - Parc Inopolis
206 route de Vouffes • 69230 Saint-Genis-Laval
Ed04-25

IMPLEMENT INNOVATIVE SOLUTIONS

OneOrtho
medical

For further information, please contact your
OneOrtho solution distributor.

OneOrtho - Parc Inopolis
206 route de Vouffes • 69230 Saint-Genis-Laval
Ed04-25

IMPLEMENT INNOVATIVE SOLUTIONS

Appendix

OneOrtho
medical

For further information, please contact your
OneOrtho solution distributor.

 OneOrtho - Parc Inopolis
206 route de Vouffes • 69230 Saint-Genis-Laval
Ed04-25

IMPLEMENT INNOVATIVE SOLUTIONS

Copyright information

This document and what it describes are protected by the law of March 11, 1957 on literary and artistic property, amended by the law of July 3, 1985, as well as by copyright laws and international conventions.

Without the written permission of ONE ORTHO no part of this manual, including the products and software described herein, may be reproduced, transcribed, stored in a database system, or translated into any language, in any form or by any means, except for documentation retained by the purchaser for backup purposes. Any person failing to comply with these provisions shall be guilty of counterfeiting and liable to the penalties provided for by law.

ONE ORTHO provides this manual "as is" without warranty of any kind, either expressed or implied, or conditions of merchantability or fitness for a particular purpose. In no event shall ONE ORTHO, its directors, officers, employees or agents be liable for any indirect, special, incidental or consequential damages even if ONE ORTHO has been advised of the possibility of such damages resulting from any defect or error in this manual or product.

Product and company names appearing in this manual may or may not be registered trademarks or subject to copyrights of their respective companies, and are used only for purposes of identification or explanation, and for the sole benefit of the owners, without any intention of infringement.

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
OneOrtho **LEFT** **Tibial varus 3.8°** Cartilage not included in CT imaging Initial HKA 176.2°

Preoperative View Planning

INFO PATIENT **Knee left** Surgical Ref: DEMO-Left Last Name: DEMO Born on: 01/01/1976 First Name: Left Surgery: 01/06/2026

Actions [Icons]

⚠️ This view is based on the acquisition of the patient's leg in the imaging system. The joint is not weight-bearing.

Angle between the femoral mechanical axis and the femoral anatomical axis in the frontal plane. Measured internally. **HKS 6.2°**

Distance between the distal end of the internal part of the condyle and the bottom of the internal tibial cup. **Internal space 1.6mm**

Angle between the femoral mechanical axis and the axis passing through the distal femoral condyles (medial and lateral) in the frontal plane. Measured externally. **LDFA 86.8°**

External space 5.5mm

MPTA (Varus) 87.2°

Angle between the tibial mechanical axis and the axis passing through the medial and lateral centers of the tibial plateau, internally, in the frontal plane.

Distance between the distal end of the external part of the condyle and the bottom of the external tibial cup. **HKA 176.2°**

Angle between the femoral mechanical axis and tibial mechanical axis in the frontal plane. Measured internally.

Angle between the tibial mechanical axis and the femoral mechanical axis in the sagittal plane (complementary to the large angle). Flexum: knees in flexion. Recurvatum: knees in extension. **Flexion 10.1°**

Sum of the lengths of the femoral mechanical axis and the tibial mechanical axis in 3D (non-projected). **833 mm**

Tibial slope 13.1°

Angle between the perpendicular to the anteroposterior median axis and the tibial mechanical axis in the sagittal plane.

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The image shows a configuration panel for a patella implant. At the top, there is a green button labeled 'PATELLA'. Below it, there are five rows of controls:

- Reference:** A dropdown menu with 'Posterior' selected.
- Meas. ref.:** A dropdown menu with 'Posterior condyles' selected.
- Range:** A dropdown menu with 'RESURFACAGE' selected.
- Diameter:** A numeric input field with '34mm' displayed, flanked by minus and plus buttons.
- Thickness:** A numeric input field with '8mm' displayed, flanked by minus and plus buttons.

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