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www.oneortho-medical.com

TABLE OF CONTENTS:

1-	GENERAL	65
2-	TERMS OF USE.....	71
3-	LIFETIME OF THE DEVICE	72
4-	GENERAL INFORMATION ON DISPLAYS	73
5-	SURGICAL PREFERENCES	77
6-	“PREOPERATIVE VIEW” PAGE.....	79
7-	PLANNING	82
7.1.	MODIFICATION OF PLANNING	82
7.1.1.	<i>Planning the femur</i>	97
7.1.2.	<i>Planning the tibia</i>	102
7.1.3.	<i>Planning the patella</i>	106
7.1.4.	<i>3D planning visualization</i>	108
7.1.5.	<i>Indication of space generated</i>	110
7.2.	PLANNING VALIDATION AND REPORT GENERATION	112
8-	MAINTENANCE.....	114
9-	CYBERSECURITY	114
10-	INCIDENT / BUG	114
11-	CONTACT	114
12-	MANUFACTURER.....	115
13-	DEVICE INFORMATION AND ASSOCIATED SYMBOLS	115

Restrictions on use:

This software is intended for use by professionals (Orthopaedic surgeons). In no way should it be used as the sole basis for clinical decisions on diagnosis, care or management of the patient. The plausibility of information obtained through the software must be systematically verified clinically before it is used to treat patients. Any application of the program's medical information that would not be in line with the original concept or intended use of the program is not recommended and will be considered as misuse of the software.

General

1- General

DEVICE NAME / TRADE NAME:

The name of the device range is "3D KNEE PLANNING SOFTWARE".

The reference presented in this document is the "Digital Planner Total" reference.

INTENDED USE:

3D KNEE PLANNING SOFTWARE is intended to provide orthopaedic surgeons with basic support for clinical decision-making through 3D visualization and planning tools for preoperative surgery planning in patients undergoing knee arthroplasty.

GENERAL DESCRIPTION:

The 3D KNEE PLANNING SOFTWARE is a planning software designed for surgical planning of knee prosthesis. It is a device that provides information to clinical management via planning tools and 3D visualization of proximal tibial, distal femoral and patellar anatomy. Mainly by that mean, it allows the planning of the positioning, size and model of total knee prosthesis in order to plan surgery alone or to plan the surgery and design a custom-made guide, if required, through the associated workflow platform. Thus, it provides the surgeon user with information that constitutes basic level of support to take a decision regarding the pre-operative planning of the surgery. It includes a 3D library with several implant models, in their different sizes, of total knee prostheses. Once the planning has been performed and validated by the surgeon user, the software provides a planning report

Using 3D items reconstructed from the DICOM images provided by the surgeon, this software displays the patient's knee (proximal tibia, distal femoral and patellar bone structure) with a 3D view and 2D views (NIfTI images).

It displays the knee prosthesis within the digitally reconstructed bone structure. The surgeon can virtually manipulate the prosthesis components (tibial, femoral, and patellar) with precision by adjusting parameters such as lateral position, femoral tibial and patellar cut depth, inclination angles, anterior notching range, and size.

The default position of the prosthesis is set by One Ortho's Digital Production Department in such a way as to reduce the case's processing time for the surgeon user and is not a planning suggestion.

SOFTWARE VERSION:

The active version of the "Digital Planner Total" software is **version 2.2**.

General

PRECAUTIONS FOR USE:

Residual risks	Recommendations	Patient consequences
The user does not understand how to open or use the 3D scene. The user is unable to find a value's definition using the "?" button which appears when hovering the value. The user cannot find the online help section. (usability risks)	Before using the device, go through the User Manual: all steps and tools are described. Also go through the Appendix 1 "Getting started Guide". Always follow good practices related to the surgery, see section "Risk control measures external to the software" stated below.	Longer operating time.
The user leaves the workstation without logging off. (usability risk) The user shares his login information. (usability risk) Access to the software by someone who should not have access rights.	Never transfer access rights to a third party, login and password are strictly personal. Always log off after you have finished using the software.	No patient consequence.
If the design of the labeling format is not carried out according to regulatory requirements (logos and information)	Make sure to understand the label. In case of doubt, consult the detailed information in the User Manual.	No patient consequence
Wrong User Manual version on the IFU website leading to software use error The user is not informed that the User Manual has been revised.	Ensure that the version of the User Manual used is identical to the version identified on the software label. Always follow good practices related to the surgery, see section "Risk control measures external to the software" stated below.	Longer operating time.

General

Residual risks	Recommendations	Patient consequences
The user uses an outdated or unverified browser. A security flaw due to using an outdated browser.	Always use a browser compliant with the manufacturer's requirements and ensure that it is not obsolete in order to avoid security failures. Always follow good practices related to the surgery, see section "Risk control measures external to the software" stated below.	Longer operating time.
Software not working. The output data (STL files) leading to the production of custom-made cutting guides does not match the planning made by the surgeon user.	Follow minimum requirements concerning hardware (including screen definition) and IT networks characteristics as defined in the User Manual to ensure correct display of device's interface as to avoid inconvenience and security issues. Always follow good practices related to the surgery, see section "Risk control measures external to the software" stated below.	No patient consequence. The medical staff must use other means to perform the plan of the surgery.
Erroneous data/Display problem	To avoid issues with the planning output data, make sure the patient's data is correct. Always follow good practices related to the surgery, see section "Risk control measures external to the software" stated below.	Error in the data. Longer operating time.
Case canceled by mistake (usability risk) Loss/forgetting of planning by the surgeon	Make sure to validate the case to generate and archive the planning report. After validation, the planning report is available on the platform.	No patient consequence. The medical staff must use other means to perform the plan of the surgery.

General

Residual risks	Recommendations	Patient consequences
Bugs, non-compliance and/or material safety issues not reported to ONE ORTHO	Please let us know if you experiment any bug (see §contact)	No patient consequence.
The user needs to use the user manual in paper form.	At the user's request, the User Manual can be supplied in paper form within 7 days of the request, at no additional cost.	No patient consequence.

Risk control measures external to the software:

The use of 3D KNEE PLANNING SOFTWARE does not replace the usual intraoperative control measures.

In order to avoid increasing operating time, the user is required to follow risk control measures external to the software:

- Verify that the data on the screen matches the patient data (name/anonymized code, date of birth, side to be operated, date of intervention...)
- When using traditional instrumentation, respect the operating technique established by the manufacturer of the prosthesis, particularly:
 - the use of traditional instrumentation dedicated to measure alignment and the needed prosthesis size.
 - the use of trial prosthesis before implanting the definitive metallic prosthesis.
- When using custom-made cutting guides, verify that the cutting planes are correct regarding the patient bony anatomy, as mentioned in the Patient-Specific Instrumentation's operating technique, by:
 - using the instrumentation to check that the cutting slots are well positioned,
 - using the alignment control rod to check the posterior tilt

The data presented by the software is not a planning suggestion, and the choice of planning is made by the surgeon user.

It is forbidden to use several web browsers or browser tabs to simultaneously access different pages of the application.

INSTALLATION SPECIFICATIONS:

Operating systems:

- Windows,
- Mac OS.

Browser:


- Google Chrome

The browser used must be compatible with WebGL technology:

<https://caniuse.com/?search=webgl>

General

The zoom level must be 100% for optimal use of the platform:

General

Screen definition:

The platform is fully responsive, i.e. it adapts to the size of the screen. This adaptation can, on small smartphone-type formats, modify the layout of the user screen. It is recommended to use a minimum of 1200px in 16:9 to display all the information in the application.

Internet connection:

An Internet connection is required to use the application. A minimum connection speed of 10 Mbps is recommended.


Equipment:

The client workstation hardware required for a comfortable user experience is:

- **Processor:** The processor must be equivalent to an Intel Core i5 (sixth-generation or higher).
- **Graphics card:** The graphics card must have a minimum video memory of 2048 MB.

ACCESSIBILITY OF DOCUMENTATION :

This User Manual can be accessed at any time by clicking on the  icon.

The regulatory label is available on the planner by hovering the mouse over the icon .

Lifetime of the device

2- Terms of use

Indications

3D KNEE PLANNING SOFTWARE is indicated to assist orthopaedic surgeons in the preoperative planning of knee arthroplasty procedures in patients with knee osteoarthritis.

Intended users

The intended users of the 3D KNEE PLANNING SOFTWARE are orthopedic surgeons with specialized knowledge in musculoskeletal anatomy, biomechanics, and reconstructive surgery, including knee arthroplasty procedures. The surgeon is responsible for reviewing and approving the preoperative surgical plan using the device.

There is no specific software training required to use the device for the first time.

Target population

3D KNEE PLANNING SOFTWARE is intended for use in male and female adult patients of advanced age who have reached bone maturity, diagnosed with knee osteoarthritis and scheduled for knee arthroplasty (age close to 65-75 years old).

The device is not intended for use in pediatric patients, pregnant and breast-feeding women, or other special populations.

Performance claims

3D KNEE PLANNING SOFTWARE provides clinical management with tools to planify knee arthroplasty surgery and predict the implant size of knee prosthesis components for individual patients undergoing knee arthroplasty. The exact femoral component size is predicted in 90% of cases, and the exact tibial component size in 75% of cases.

Clinical benefits

The clinical benefits of 3D KNEE PLANNING SOFTWARE are indirect. By enabling accurate preoperative planning, the device supports surgical workflow efficiency. Its use has been associated with a reduction in total procedure time by approximately 15% compared to conventional knee arthroplasty

Contraindications

The planning software must not be used:

- For joints other than the knee joint.
- For planning other prostheses ranges than those referenced in the library.
- When knee arthroplasty is contraindicated.

Adverse side effects

No adverse side effects have been recorded in association with the device.

Lifetime of the device

3- Lifetime of the device

This software will be maintained for two years.


During this period, One Ortho guarantees the performance of the software and undertakes to provide updates to enable its optimal use.

General information on displays



4- General information on displays

In this section, we'll describe all the information shown and the possible interactions in this insert, which is common to both pages of the application. As required, the user can navigate from the preoperative page to the planning page and vice-versa.


Whichever page is displayed, the application's top banner looks like this:

To move from one page to another, simply click on the corresponding navigation arrows:


-  The "Preoperative view" page is active.
-  The "Planning" page is active.

In the top right-hand corner of the application, the following information is displayed:


General information on displays


Below the navigation arrows, the following information is displayed:

Tools for specific actions are also available (can only be used on the Planning page):


Pressing the Colourimetry button  will display two sliders that can be used to adjust the contrast and brightness:


Pressing the button  will take a screenshot of part of the interface, which can be saved to the user's computer or added to the planning report. To do this, select a rectangular area by holding down the left mouse button and moving the mouse, while remaining within the browser window.

General information on displays


Once the button is released, the following interface is displayed to name the screenshot. Please note that if the selected area extends beyond the browser, the screenshot will be blank.

The buttons at the bottom allow you to cancel, save the image to the user's workstation, or add it to the planning report.

Pressing the button  will display the following interface, allowing you to choose the colours of the bones and implants.


Choice of bone colour

Choice of implant colour
from a range of predefined
colours


Free choice of implant
colour

General information on displays

Tools for managing navigation in the application are available (only the temporary save button is active on the "Preoperative view" page):

Finally, tools for managing displays, accessing explanatory documentation and disconnecting from the platform are available at any time during the planning process:


Surgical preferences

5- Surgical preferences

Users can specify their preferred planning parameters.


To do this, access the preference settings screen via the  button:

The interface for this page is as follows:


The user can then set the following parameters:

- The appearance of the planning theme:


- The default page displayed when the case is opened

Surgical preferences

- **Default values for the display of cut bones (femur, tibia, and patella):**
 - Hidden
 - Displayed
 - If displayed:
 - Transparency level (0% = transparent, 100% = opaque)
 - Option to display cut surfaces in color
- **Default values for the display of uncut bones (femur, tibia, and patella):**
 - Hidden
 - Displayed
 - If displayed:
 - Transparency level (0% = transparent, 100% = opaque)
- **Default values for certain femoral planning information:**
 - Reference
 - Rotation axis
 - Femoral implant range
- **Default values for the display of the femoral implant:**
 - Hidden
 - Displayed
 - If displayed:
 - Transparency level (0% = transparent, 100% = opaque)
- **Default values for certain tibial planning information:**
 - Reference
 - Tibial implant range
- **Default values for the display of the tibial implant:**
 - Hidden
 - Displayed
 - If displayed:
 - Transparency level (0% = transparent, 100% = opaque)


To validate any changes made and return to the planning, the user must click on the button

Preoperative view page


6- "Preoperative view" page

One of the application's pages is entitled "Preoperative view". It contains all data relating to the patient's native anatomy:

The views shown are illustrations based on the DICOM images of the patient's leg according to the established imaging protocol. The joint is not under load.


Preoperative view page

Definition of displayed measurements:

- **HKS**: Angle between the femoral mechanical axis and the femoral anatomical axis in the frontal plane. Measured internally.
- **HKA**: Angle between the femoral mechanical axis and tibial mechanical axis in the frontal plane. Measured internally.
- **LDFA**: Angle between the femoral mechanical axis and the axis passing through the distal femoral condyles (medial and lateral) in the frontal plane. Measured externally.
- **MPTA (Tibial varus)**: Angle between the tibial mechanical axis and the axis passing through the medial and lateral centers of the tibial plateau, internally, in the frontal plane.
- **INTERNAL SPACE**: Distance between the distal end of the internal part of the condyle and the bottom of the internal tibial cup.
- **EXTERNAL SPACE**: Distance between the distal end of the external part of the condyle and the bottom of the external tibial cup.

Preoperative view page



Definition of displayed measurements:

- **Flexion:** Angle between the tibial mechanical axis and the femoral mechanical axis in the sagittal plane. If the value is positive, the knee is in flexion. If it is negative, the knee is in extension.
- **Tibial slope:** Angle between the perpendicular to the anteroposterior median axis and the tibial mechanical axis in the sagittal plane.

Preoperative view page

At any time, the user can view the measurement definition by clicking on it.

They also may view all the definitions by pressing the  button:

7- Planning

7.1. Modification of planning


One of the pages in the application is entitled "Planning". It contains all the tools needed to plan the current case:

Planning

Dynamic 3D visualization
and 2D visualization
layout insert

Interactive insert for 2D visualization,
planning, movement and positioning
of the femoral implant

Visualization tools
access panel

Dynamic insert to indicate generated spaces

Dynamic HKA angle
insert calculated in real
time

Interactive insert for 2D
visualization, planning,
movement and positioning of the
tibial implant


For patella planning , you must first activate the patella visualization (by double clicking)
The screen then becomes this one:

Interactive visualization
insert for patellar
planning

Interactive patellar
implant planning insert

Button to press by
double clicking to access
to the patellar tools

Planning


Description of the visualization tools access panel:

Scanner view/3D view

Scanner view
active

“3D model visualization” mode activated:

Planning

“Scanner image visualization” mode activated:

In this scanner view, the user can move section by section using the mouse wheel and view the overlaying of the selected sections in 3D mode and the selected implants. In this way, the user can see the coverage of the selected implant and validate or not his choice of planning.

A zoom slider is available for each slice in the CT scan view. The user can drag the slider upward to zoom in or downward to zoom out. The zoom is applied in real time to the corresponding slice, with no view reload required.

This feature allows the user to inspect anatomical details more closely, particularly when assessing bone structures during tibial stem extension planning.

When a slice is zoomed in, the area of interest may fall outside the visible frame. The user can navigate within the enlarged image without changing the zoom level by holding down the right mouse button on the slice and moving the mouse in the desired direction. The slice follows the mouse movement. The user releases the right mouse button to stop panning.

The zoom level remains unchanged throughout the panning operation.

Planning

Display of the Patella and planning

Display
Patella
vizualization
tools


Planning of the
patella


Patella hidden

Display Patella vizualization tools:


Planning of the patella:

Patella hidden:

Planning

Visualization "Bone cut models"

**Bone cut
models
displayed**


**Bone cut models
displayed
transparently**


**Bone cut models
hidden**

Bone cut models displayed :

Bone cut models displayed transparently:

Planning


Bone cut models hidden:


Planning


Visualization "Femoral implant"

Femoral implant displayed:


Femoral implant displayed transparently:

Femoral implant hidden:

Planning

Vizualization « Patellar implant »

Patellar
implant
displayed


Patellar implant
displayed
transparently


Patellar implant

Patellar implant displayed :


Patellar implant displayed transparently :

Planning

Patellar implant hidden :

Visualization “Tibial implant”


**Tibial implant
displayed**

**Tibial implant
displayed
transparently**


**Tibial implant
hidden**

Tibial implant displayed:

Planning


Tibial implant displayed transparently:

Planning

Tibial implant hidden :

Visualization "Bone resections"

**Bone
resections
displayed**


**Bone resections
displayed
transparently**


**Bone resections
hidden**

Bone resections displayed:


Bone resections displayed transparently:

Planning

Bone resections hidden:


Planning

Change in sagittal view orientation


Left case:
Lateral view
Right case:
Medial view

Left case: Medial view
Right case: Lateral view


Medial sagittal view (left case):


Lateral sagittal view (left case):

Displaying/hiding axes


Axes displayed


Axes hidden

Axes displayed:

Planning


Axes hidden:

Enabling/Disabling coloured cuts


Colourful cuts disabled


Colourful cuts displayed

Colourful cuts displayed :

Planning

Colourful cuts disabled :

7.1.1. Planning the femur

The user is provided with a range of tools for planning the femur. They are identified on the planner by the "Purple" color. The tools for positioning and sizing the femoral implant are as follows:

Planning

FEMUR


Reference	Posterior
Rot. axis	Diaphyseal axis
Meas. ref.	Posterior condyles
Range	CONDYLE PS
Size	- 4 +

The user can then select the following parameters:

<u>Reference</u>	The choice of this reference will define the default anteroposterior (AP) positioning of the condyle.
<u>Rotation axis</u>	The choice of the axis of external rotation will determine how the condyle rotates in the horizontal plane.
<u>Measurement reference</u>	The choice of this reference will define the 0 value for external rotation according to the chosen reference frame.
<u>Femoral implant range</u>	One of the distributor's ranges
<u>Size</u>	Guaranteed size in the ranges offered by the distributor


Planning

The following tools are available for making femoral cuts:

The user can then modify all these parameters to plan the femoral cuts for their case.

The indicative measurements displayed induced by the cut settings are as follows:


Planning


The display of femoral views can be modified according to the user's wishes. To do so, click on

to the right of the planner:

This will open a hidden window containing the 3D planning view and the settings for the various views:

Visualization of the distal femoral view

Planning


Visualization of the frontal femoral view

Visualization of the sagittal femoral view


These views can be combined by clicking one after the other.

For example, the user can choose to view only the distal and frontal femoral views:


The menus for sizing/positioning the femoral implant and creating femoral sections are displayed with the views selected by the user. This view is closed by clicking again on the button.

In distal femoral view only, the user can position the implant on the bone model. Passing the mouse over the femoral implant in distal view will change its appearance:

By clicking on the implant in this view, the user can move it medio-laterally. The implant then takes on a new appearance during this move:

Planning

This displacement is visible in all femoral views.

7.1.2. Planning the tibia

The user is provided with a range of tools for planning the tibia. They are identified on the planner by the "Blue" color. Tools for positioning and sizing the tibial implant include:

TIBIA


Reference	Posterior plateau ed... ▾
Meas. ref.	Posterior plateau ed... ▾
Range	EMBASE FIXE ▾
Size	- 6 +
Insert	6/7

The user can then select the following parameters:

<u>Reference</u>	The choice of this reference will define the default positioning of the tibial component in rotation on the tibial section.
<u>Measurement reference</u>	The choice of this reference will define the 0 value of the tibial component rotation according to the chosen reference.
<u>Tibial implant range</u>	One of the distributor's ranges
<u>Size</u>	Guaranteed size within the range offered by the distributor (+/- one size of that chosen for the condyle)
<u>Insert</u>	This is an indication of the plastic spacer size compatible with the selected size.


The tools for making tibial cuts are as follows:

Planning

The user can then modify all these parameters to plan the tibial cuts for their case.


The indicative measurements displayed induced by the cut settings are as follows:


Planning

The display of tibial views can be modified according to the user's wishes. To do so, to the

right of the planner, click on >> :

This will open a hidden window containing the 3D planning view and the settings for the various views:

Planning

Visualization of the proximal tibial view


Visualization of the frontal tibial view

Visualization of the sagittal tibial view

These views can be combined by clicking one after the other.


For example, the user can choose to view only the proximal and frontal tibial views:


The menus for tibial implant sizing/positioning and tibial sectioning are displayed according to the views selected by the user. This set view is closed by clicking again on the button .

In the distal tibial view only, the user can position the implant on the bone model. Passing the mouse over the tibial implant in distal view will change its appearance:

Planning

By clicking on the implant in this view, the user can move it medio-laterally. The implant then takes on a new appearance during this move:

This displacement is visible in all tibial views.

7.1.3. Planning the patella

The user has in his possession various tools in order to carry out the planning of the patella. They are identified on the planner by the color "Green".

The tools related to the positioning and sizing of the patellar implant are as follows:

 PATELLA

Reference

Posterior

Meas. ref.

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Range

RESURFACAGE

Diameter

-

34mm

+

Thickness

-

8mm

+

The user will then be able to choose the following parameters:

<u>Reference</u>	The choice of this reference will define the default anterior/posterior positioning of the patellar implant.
<u>Measurement Reference</u>	The choice of this reference will set the value 0 of the external rotation according to the chosen reference. This reference is common to that of the femur.
<u>Patellar implant range</u>	One of the distributor's ranges

Planning

Diameter


The choice of the diameter of the patellar implant from the available sizes offered by the distributor.

Thickness


The choice of patellar implant thickness from the available sizes offered by the distributor.

The tools for making patellar cuts are as follows:


Rotation of the patellar implant in the axial plane

Rotation of the patellar implant in the sagittal plane

Translation of the patellar implant in proximal/distal plane


The user will then be able to modify all these parameters to plan the patellar resections of his case.

The indicative measurements displayed induced by the adjustment of the cuts are as follows:

Patellar offset : distance between the prosthetic patellar apex and the prosthetic trochlear floor

Thickness of the patella after resection and thickness of the resection


Remaining patellar thickness: in extension, the distance between patellar resection and femoral resection

Planning


7.1.4. 3D planning visualization

Throughout the planning process, the user can view the implant sizing and position on the 3D bone models of the current case. To do so, click on the  button on the right-hand side of the planner:


Planning

The following 3D view is displayed:


Users can move the entire 3D model by holding down the left mouse button on the scene presenting the femur, tibia and patella.

They can also use the following visualization tools:

Visualization of bone models in frontal view

Visualization of bone models in sagittal view

Planning

Display/Hide the femur

Display/Hide the tibia

Display/Hide the patella


7.1.5. Indication of space generated

On the planning page, the user has access to the space values generated by the cuts:


The sum of the femoral and tibial sections gives the space generated by the sections, and thus the implant thickness to be added. Each measurement is the sum of the internal or external cuts, in two different configurations: flexion or extension. The four sums obtained are as follows (illustrated below):

- **Internal space in flexion:** Internal posterior femoral section + Internal proximal tibial section
- **External space in flexion:** External posterior femoral section+ External proximal tibial section
- **Internal space in extension:** distal medial femoral section + proximal medial tibial section
- **External space in extension:** distal external femoral section+ proximal external tibial section

Planning

Planning


7.2. Planning validation and report generation

Once the surgeon has set up all their cuts, dimensioned and positioned the implants, the user must validate definitively his planning. As a result, the user will no longer be able to modify what he been done, and the validated output data will be used by ONE ORTHO to produce the associated custom guides.

To do this, the surgeon clicks on the button in the upper menu of the planning screen. They must then confirm the validation of their case by clicking on "Yes" in the window that appears:

Planning

Once validated, a report summarizing the planning carried out will be issued to the surgeon.
The user is redirected to the procedure follow-up.

Maintenance, Cybersecurity & Contact

8- Maintenance

A yearly update is carried out for safety maintenance, automatically and completely transparently to the user.

9- Cybersecurity

Our approach to managing cybersecurity is based on risk management at every stage of the software lifecycle, integrating cybersecurity concepts right from the design stage. Security measures are applied to user authentication to the web application, which is managed by a login and password that gives only one user access to his or her own cases, thus restricting access to patient data to authorized persons only. We therefore ask you to be careful not to divulge your password, and to log out when you stop using the application.

All cybersecurity-related risks are integrated and controlled in the medical device's risk management file.

10- Incident / bug

Despite the rigor with which we have developed this software, any serious incident involving the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

A list of known bugs is available at the following link :

<https://github.com/oneorthomedical/3D-KNEE-PLANNING>

Should a bug not listed here be detected, please forward the information to ONE ORTHO.

11- Contact

Contact :

If you have any questions about how to use the platform, please contact ONE ORTHO at contact@oneortho-medical.com or +33 4 26 78 76 74.

The ONE ORTHO team is available to answer your questions from Monday to Friday, 9am to 6pm.

Maintenance, Cybersecurity & Contact

12- Manufacturer

ONE ORTHO
Parc INOPOLIS,
206 Route de Vourles
69230 SAINT-GENIS-LAVAL
FRANCE
contact@oneortho-medical.com

Class IIa medical device, manufactured by ONE ORTHO, compliant with the requirements of EU Regulation 2017/745.

Year of 1st CE marking (under EU Regulation 2017/745): 2026
2025 One Ortho. All rights reserved.

13- Device information and associated symbols

Light theme:

 Manufacturer	 Date of manufacture	 Unique device identifier	 Medical device
 Consult the User Manual	 Software version		

Dark theme:

 Manufacturer	 Date of manufacture	 Unique device identifier	 Medical device
 Consult the User Manual	 Software version		

Maintenance, Cybersecurity & Contact

Copyright information

This document and what it describes are protected by the law of March 11, 1957 on literary and artistic property, amended by the law of July 3, 1985, as well as by copyright laws and international conventions.

Without the written permission of ONE ORTHO no part of this manual, including the products and software described herein, may be reproduced, transcribed, stored in a database system, or translated into any language, in any form or by any means, except for documentation retained by the purchaser for backup purposes. Any person failing to comply with these provisions shall be guilty of counterfeiting and liable to the penalties provided for by law.

ONE ORTHO provides this manual "as is" without warranty of any kind, either expressed or implied, or conditions of merchantability or fitness for a particular purpose. In no event shall ONE ORTHO, its directors, officers, employees or agents be liable for any indirect, special, incidental or consequential damages even if ONE ORTHO has been advised of the possibility of such damages resulting from any defect or error in this manual or product.

Product and company names appearing in this manual may or may not be registered trademarks or subject to copyrights of their respective companies, and are used only for purposes of identification or explanation, and for the sole benefit of the owners, without any intention of infringement.

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